



# Strips 12 tappi piani anti-impronta extra robusti


Cf. 80 strips da 12 tappi - Strips di 12 tappi ottici piani anti-impronta



## Product codes:

Reference: 12 864 4011

## Product gallery:



## Product description:

Strips di 12 tappi ottici piani anti-impronta

La superficie rientrante evita che l'area ottica venga "contaminata" da impronte accidentali

Utilizzabili in tutte le procedure in PCR

Pareti sottili ed uniformi, ottimi per chiudere le piastre e le strips per le quali sono indicate come accessorio.

Ottimizzate per PCR

Cf. 80 strips da 12 tappi

## Product features:

FLM - TYPE: the strip of caps

FLM - COMPATIBILITY: strip PCR

FLM - COMPATIBILITY: plates PCR



FLM - COMPATIBILITY: plates Real-Time PCR  
FLM - COMPATIBILITY: strip PCR Real-Time  
FLM - MATERIAL: PP  
FLM - WELLS: 12  
FLM - PURITY: DNA / RNA / DNase / Rnase / Pyrogen FREE  
FLM - STERILITY: no  
FLM - SHAPE: reentrant  
FLM - SHAPE: strip sturdy caps  
FLM - TEMPERATURE MAX: 100°C  
FLM - MINIMUM TEMPERATURE: -30°C  
FLM - OPTICAL AREA: 12.6mm  
FLM - AREA WRITABLE: you  
FLM - AUTOCLAVABLE: you  
FLM - CODING COLOR: no  
FLM - LABEL: no  
FLM - COMPLIANCE: not subject to CE  
FLM - FEATURES: strips of 12 caps-profile re-entrant anti-interference, wide-area optical type extra-robust  
FLM - SUITABLE FOR: molecular biology  
FLM - SUITABLE FOR: PCR  
FLM - SUITABLE FOR: PCR Real Time  
FLM - SUITABLE FOR: PCR Real Time FAST  
FLM - SUITABLE FOR: storage  
FLM - SUITABLE FOR: storage at -20°C  
FLM - COMPATIBLE WITH: Applied Biosystems-9800 Fast Thermal Cyclers-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-SimpliAmp™-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cyclers-96 x 0.1ml  
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Thermal Cyclers-96 x 0.2ml  
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cyclers-96 x 0.1ml  
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cyclers-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml



FLM - COMPATIBLE WITH: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast  
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml  
FLM - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-Opticon 2 System-96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-Opticon System-96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml  
FLM - COMPATIBLE WITH: Stratagene-Mx3000P™-96 x 0.2ml  
FLM - COMPATIBLE WITH: Stratagene-Mx3005P™-96 x 0.2ml  
FLM - COMPATIBLE WITH: Stratagene-Mx4000™-96 x 0.2ml  
FLM - COMPATIBLE WITH: Stratagene-Robocycler Gradient 96-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TOptical-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml  
FLM - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml