



Strips 12 pozzetti qPCR da 0,2ml

Cf. 800 pezzi - Strips 12pozzetti da 0,2ml per qPCR, profilo regolare

Product codes:

Reference: 12 954 7841

Product gallery:



Product description:

Strips da 0,2ml, profilo regolare, eccellenti per PCR Real-Time

La scelta migliore per gli apparecchi ABI e BioRad (non FAST) per PCR Real-Time che accettano provettine da 0,2ml a profilo regolare.

Non compatibile con apparecchi FAST.

Cf. 800pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Leggenda:

Grado qPCR 1: sviluppato per segnale ottimo in PCR Real-Time

Grado qPCR 2: sviluppato per avere il segnale al top, esente da interferenze, in PCR Real-Time



Product features:

SP - TYPE: strip 12 microtubes
SP - CAP: senza
SP - VOLUME: 0.2ml
SP - ABILITY: 0.2ml
SP - WELLS: 12
SP - PROFILE: adjust
SP - HEIGHT: 20.9 / 21.8mm ca
SP - SHAPE: strip sturdy
SP - FUND: tapered
SP - MATERIAL: PP
SP - APPLICATIONS: PCR
SP - APPLICATIONS: PCR Real Time
SP - AUTOCLAVABLE: you
SP - COLOR: neutral
SP - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITY: no
SP - CONFORMITY: not subject to CE
SP - MAX TEMPERATURE: 100°
SP - MINIMUM TEMPERATURE: -30°C
THE SP - PACKAGE: envelope
SP - FEATURES: strip 12x0.2 for Real Time PCR and PCR, extra robust
SP - OPTIONAL ACCESSORIES: the strip of caps
SP - OPTIONAL ACCESSORIES: the grid assembly blue
SP - SUITABLE FOR: PCR
SP - SUITABLE FOR: PCR Real Time
SP - SUITABLE FOR: storage
SP - WALLS: thin 0.30 mm
SP - TREATMENT SIGNAL: walls of frosted
SP - SIGNAL-STRENGTH: 125
SP - POSSIBLE INTERFERENCE SIGNAL: 70
SP - RESISTANT: freezing at -20°C
SP - EVAPORATION SAMPLE: the lesser of 3%
SP - BINDING: DNA minor 0.3%
SP - BINDING: proteins less 0.8%
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well



SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-310-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3100-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3130-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3500 Genetic Analyzer-8-capillary
SP - COMPATIBLE WITH: Applied Biosystems-3700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3730-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
SP - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
SP - COMPATIBLE WITH: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml



SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-SimpliAmp™-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-Solid3-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
SP - COMPATIBLE WITH: Stratagene-Mx3000P™-96 x 0.2ml
SP - COMPATIBLE WITH: Stratagene-Mx3005P™-96 x 0.2ml
SP - COMPATIBLE WITH: Stratagene-Robocycler Gradient 96-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TOptical-96 x 0.2ml