



Strips 12 tappi ottici piani anti-impronta per Real Time PCR


Strips di 12 tappi ottici piani anti-impronta specifici per PCR Real-Time



Product codes:

Reference: 12 864 1951

Product gallery:



Product description:

Strips di 12 tappi ottici piani anti-impronta specifici per applicazioni in PCR Real-Time

La superficie rientrante evita che l'area ottica venga "contaminata" da impronte accidentali

Utilizzabili in tutte le procedure in PCR Real Time

Pareti sottili ed uniformi, ottimi per chiudere le piastre e le strips per le quali sono indicate come accessorio.

Ottimizzate per qPCR

Cf. 80 strips da 12 tappi

**Product features:**

FLM - TYPE: the strip of caps
FLM - COMPATIBILITY: strip PCR
FLM - COMPATIBILITY: plates PCR
FLM - COMPATIBILITY: plates Real-Time PCR
FLM - COMPATIBILITY: strip PCR Real-Time
FLM - MATERIAL: PP
FLM - WELLS: 12
FLM - PURITY: DNA / RNA / DNase / Rnase / Pyrogen FREE
FLM - STERILITY: no
FLM - SHAPE: reentrant
FLM - SHAPE: strip sturdy caps
FLM - TEMPERATURE MAX: 100°C
FLM - MINIMUM TEMPERATURE: -30°C
FLM - OPTICAL AREA: 12.6mm
FLM - AREA WRITABLE: you
FLM - AUTOCLAVABLE: you
FLM - CODING COLOR: no
FLM - LABEL: no
FLM - COMPLIANCE: not subject to CE
FLM - FEATURES: strips of 12 caps-profile re-entrant anti-interference, wide-area optical type robust
FLM - SUITABLE FOR: molecular biology
FLM - SUITABLE FOR: PCR
FLM - SUITABLE FOR: PCR Real Time
FLM - SUITABLE FOR: PCR Real Time FAST
FLM - SUITABLE FOR: storage
FLM - SUITABLE FOR: storage at -20°C
FLM - COMPATIBLE WITH: Applied Biosystems-310-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-3100-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-3130-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-3700-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-3730-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBLE WITH: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBLE WITH: Applied Biosystems-9800 Fast Thermal Cyclers-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast



FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
FLM - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-SimpliAmp™-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-Solid3-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
FLM - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml



FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
FLM - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-Opticon System-96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBLE WITH: Roche-cobas z 480-96 x 0.2ml
FLM - COMPATIBLE WITH: Roche-LightCycler ® 480 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE WITH: Stratagene-Mx3000P™-96 x 0.2ml
FLM - COMPATIBLE WITH: Stratagene-Mx3005P™-96 x 0.2ml
FLM - COMPATIBLE WITH: Stratagene-Mx4000™-96 x 0.2ml
FLM - COMPATIBLE WITH: Stratagene-Robocycler Gradient 96-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TOptical-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
FLM - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml