




Piastra sub-skirted, basso profilo, specifiche ABI FAST / StepOnePlus


Cf. 25 pezzi - Piastre sub-skirted, basso profilo, specifiche per ABI Fast e StepOne Plus

Product codes:

Reference: 12 893 1981



Product gallery:



Product description:

Piastre da 0,1ml, basso profilo, eccellenti per PCR Real-Time. Colore neutro

Compatibilità limitata e specifica agli apparecchi ABI Fast e StepOne Plus. Facilmente tagliabili.

Cf. 25pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Leggenda:

Grado qPCR 1: sviluppato per segnale ottimo in PCR Real-Time

Grado qPCR 2: sviluppato per avere il segnale al top, esente da interferenze, in PCR Real-Time



Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.1ml
PSC - FORM: sub-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V
PSC - CLOSING: through film
PSC - CLOSING: via caps
PSC - CLOSING: mat
PSC - GRID: alphanumeric
PSC - APPLICATIONS: PCR Real-Time FAST
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: neutral
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - STACKABLE: you
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - MAIN FEATURES: plate PCR 96-well 0.1-PCR Real-Time FAST skirted specific ABI
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - SUITABLE FOR: PCR Real-Time FAST
PSC - SUITABLE FOR: sample preparation
PSC - SUITABLE FOR: storage
PSC - PROFILE: low type FAST/0.1
PSC - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
PSC - CAPACITY: 0.2ml
PSC - PACKAGE: bulk
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%



Shop name

PSC - BINDING: proteins less 0.8%