



Shop name

Microtubi PCR da 0,2ml tappo a cupola anti-contaminazione



Cf. 1000 pezzi - Microtubi da 0,2ml per PCR, tappo a cupola, anti-contaminazione

Product codes:

Reference: 12 879 9511



Product gallery:



Product description:

Microtubi da 0,2ml, tappo a cupola, anti-contaminazione

Compatibili con molti apparecchi per PCR che accettano provettine da 0,2ml con tappo a cupola.

Profilo regolare.

Cf. 1000pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Product features:

MTB - CAP: the semi-dome



MTB - CAP: the semi-dome of contamination
MTB - COMPLIANCE: not subject to CE
MTB - TYPE: PCR
MTB - VOLUME: 0.2ml
MTB - CAPACITY: 0.2ml
MTB - CAP FORABILE: no
MTB - INSERTED CAP: no
MTB - OPTICAL AREA: 6.8mm
MTB - FORM: adjust
MTB - FUND: tapered
MTB - SELF STANDING: no
MTB - MATERIAL: PP
MTB - SIZED: no
MTB - OPACITY: ultra-transparent
MTB - AREA WRITABLE: no
MTB - TEMPERATURE RANGE: -30°C a +100°C
MTB - TEMPERATURE MAX: 100°C
MTB - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - AUTOCLAVABLE: you
MTB - STERILITY: no
MTB - SIZE: h 20.9/21.8mm
MTB - FEATURES: microtubes 0.2 ml with cap dome anti-contamination
MTB - SUITABLE FOR: sample preparation
MTB - SUITABLE FOR: storage
MTB - SUITABLE FOR: PCR
MTB - SUITABLE FOR: PCR Real Time
MTB - SUITABLE FOR: PCR Real Time FAST
MTB - RESISTANT: freezing at -20°C
MTB - PACKAGING: envelope bulk
MTB - WALLS: thin 0.3 mm
MTB - EVAPORATION SAMPLE: the lesser of 3%
MTB - BINDING: DNA minor 0.3%
MTB - BINDING: proteins less 0.8%
MTB - SIGNAL-STRENGTH: 100
MTB - INTERFERENCE SIGNAL: 100
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml



MTB - COMPATIBLE WITH: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml



MTB - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml