



Microtubi PCR da 0,5ml tappo piatto, graduati

Cf. 1000 pezzi - Microtubi da 0,5ml per PCR, tappo piatto, graduati

Product codes:

Reference: 12 879 3051



Product gallery:



Product description:

Microtubi da 0,5ml, tappo piatto, graduati.

Compatibili con molti apparecchi per PCR che accettano provettine da 0,5ml con tappo piatto.

Profilo regolare.

Cf. 1000pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Product features:

MTB - CAP: dish

MTB - COMPLIANCE: not subject to CE

MTB - TYPE: PCR



MTB - TYPE: Eppendorf
MTB - VOLUME: 0.5ml
MTB - CAPACITY: 0.5ml
MTB - CAP FORABILE: no
MTB - INSERTED CAP: no
MTB - OPTICAL AREA: 5mm
MTB - FORM: adjust
MTB - FUND: tapered
MTB - SELF STANDING: no
MTB - MATERIAL: PP
MTB - SIZED: you
MTB - OPACITY: ultra-transparent
MTB - AREA WRITABLE: no
MTB - TEMPERATURE RANGE: -30°C a +100°C
MTB - TEMPERATURE MAX: 100°C
MTB - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - AUTOCLAVABLE: you
MTB - STERILITY: no
MTB - SIZE: h 30mm ca
MTB - FEATURES: microtubi da 0.5ml con tappo piatto graduati per PCR e PCR Real-Time
MTB - SUITABLE FOR: sample preparation
MTB - SUITABLE FOR: storage
MTB - SUITABLE FOR: PCR
MTB - SUITABLE FOR: PCR Real Time
MTB - RESISTANT: freezing at -20°C
MTB - PACKAGING: envelope bulk
MTB - WALLS: thin 0.3 mm
MTB - BINDING: proteins less 0.8%
MTB - BINDING: DNA minor 0.3%
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-48 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 77 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-20 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-20 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TGradient-48 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TPersonal-20 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional-60 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard-60 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-30 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-18 18 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TRIO Thermoblock-20 x 0.5ml



MTB - COMPATIBLE WITH: Whatman Biometra-UNO II-48 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-40 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-GeneCycler-12 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 480-48 x 0.5ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-60 x 0.5ml
MTB - COMPATIBLE WITH: Applied Biosystems-Veriti™ 60-Well Thermal Cycler-60 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-60 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 12 x 0.5ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-60 x 0.5ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 77 x 0.5ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 77 x 0.5ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler Personal-25-16 16 x 0.5ml
MTB - COMPATIBLE WITH: Eppendorf-Thermomixer C-24 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-100-60 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-150 MiniCycler-16 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-60 Alpha unit 60 x 0.5ml