



Strips 8 microtubi da 0,2ml modello Classic

Cf. 120 strips - Strips di 8 microtubi da 0,2ml per PCR modello "Classic" profilo regolare

Product codes:

Reference: 12 954 7781



Product gallery:



Product description:

Strips 8 microtubi da 0,2ml

Testati lotto per lotto Nuclease, Endotossine e DNA Free

Compatibili con:

- Applied Biosystems: 2400, 2700, 2720, 5700, 7000, 7300, 7500, 7700, 9600, 9700, Veriti
- BioRad: mod. GeneCycler, IcyCler, IQ5, MyCycler, MyIQ, Ptc-100, Ptc-200, Ptc-2250, MiniOpticon, Chromo4, CFX96
- Eppendorf: mod. Mastercycler, MC Gradient, MC ep Gradient, MC ep Gradient S, MC Pro, MC ep Realplex
- MJ Research: mod. PTC-100, PTC-200, PTC-225
- MWG: mod. Primus 96



-
- Operon: mod. OPR-PCR96
 - Stratagene: mod. Robocycler, MX3000P, MX3005P, MX4000
 - Takara: mod. TP-3000
 - Techne: mod. Touchgene, Flexigene, Genius, Genius Quad
 - Thermo Hybaid: mod. PCR Sprint, PCR Express, Ommigene, OMN-E, Touchdown
 - Thermolyne: mod. Amplitron
 - Whatman Biometra: mod. Tpersonal, T1, Tgradient, T3000, Tprofessional, Trobot
 - Quanta Biotech: mod. AutoQ, QB96
 - Roche: NO
 - Abbot: mod. M2000rt
 - Analytik-Jena: mod. FlexCycler
 - Bioer: Life Express, Little Genius, XPcycler
 - Corbett Research: PalmCycler
 - Esco: mod. Swift Maxi
 - Dupont: mod. BAQ Q7
 - EuroClone: mod. One**, peqSTAR
 - Finnzymes mod Piko -24
 - Gene Technologies mod. G-Storm
 - Labnet mod. MultiGene
 - LongGene mod. MyGene
 - Peqlab mod. PeqSTAR, Primus 96
 - Scinics mod. EZ Cycler
 - Sensoquest mod. Labcycler



- TaKaRa mod. PCR Thermal Cycler Dice

Product features:

SP - TYPE: strip da 8 microtubi
SP - CAP: senza
SP - VOLUME: 0.2ml
SP - ABILITY: 0.2ml
SP - WELLS: 8
SP - PROFILE: adjust
SP - HEIGHT: 20.9 / 21.8mm ca
SP - SHAPE: striscia sottile
SP - FUND: tapered
SP - MATERIAL: PP
SP - APPLICATIONS: PCR
SP - APPLICATIONS: PCR Real Time
SP - AUTOCLAVABLE: you
SP - COLOR: neutral
SP - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITY: no
SP - CONFORMITY: not subject to CE
SP - MAX TEMPERATURE: 100°
SP - MINIMUM TEMPERATURE: -30°C
THE SP - PACKAGE: envelope
SP - FEATURES: strip da 8x0.2 per PCR Real Time e PCR ultra-trasparenti
SP - SUITABLE FOR: PCR
SP - SUITABLE FOR: PCR Real Time
SP - SUITABLE FOR: storage
SP - WALLS: thin 0.30 mm
SP - TREATMENT SIGNAL: nessuno
SP - SIGNAL-STRENGTH: 100
SP - POSSIBLE INTERFERENCE SIGNAL: 100
SP - RESISTANT: freezing at -20°C
SP - EVAPORATION SAMPLE: the lesser of 3%
SP - BINDING: DNA minor 0.3%
SP - BINDING: proteins less 0.8%
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml



SP - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml



SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
SP - COMPATIBLE WITH: Stratagene-Mx3000P™-96 x 0.2ml
SP - COMPATIBLE WITH: Stratagene-Mx3005P™-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 30 x 0.5ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 30 x 0.5ml