



copy of Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted

Cf 25pz - Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted. Cf 25pz



Product codes:

Reference: 12 893 0701

Product gallery:



Product description:

Piastra da 0.2ml semi skirted. Utilizzabile su sistemi PCR e PCR Real Time a profilo regolare. Questa piastra è inoltre perfettamente compatibile con sistemi ABI e Life Technologies

- Colore bianco
 - Regular Profile
 - Sub skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.1ml
PSC - FORM: skirted
PSC - OPENING: normal
PSC - WELLS: 96



PSC - BOTTOM: a V
PSC - CLOSING: through film
PSC - CLOSING: via caps
PSC - CLOSING: mat
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: transparent
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - PACKAGE: bulk
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: adjust
PSC - COMPATIBLE WITH: Agilent Technologies AriaDx Realtime PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Agilent Technologies AriaMx Realtime PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena Alpha SC 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler Monoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TOne 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TOne 96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Axygen® MaxyGene™ II 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GenePro B-96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer LifePro TC-96/G/H(b)A 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Line-Gene 9600 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Opus Deepwell 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBLE WITH: Blue-Ray Biotech TurboCycler TCST-9610 96 x 0.2ml
PSC - COMPATIBLE WITH: Blue-Ray Biotech TurboCycler TCST-9620 96 x 0.2ml



PSC - COMPATIBLE WITH: Bulldog Bio GenePro™ Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GeneTouch™ Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio LifePro™ Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio LifeTouch™ Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Corbett Research Palm-Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME GenePro Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME LifeExpress B-96D 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME LifeExpress TC-96 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50a 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50i 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50l 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50s 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50t 384 Microplate
PSC - COMPATIBLE WITH: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ Maxi 96 x 0.2ml
PSC - COMPATIBLE WITH: gradient
PSC - COMPATIBLE WITH: Esco Swift™ MaxPro 96 x 0.2ml
PSC - COMPATIBLE WITH: GE Healthcare MegaBACE 1000 96 x 0.1ml
PSC - COMPATIBLE WITH: GE Healthcare MegaBACE 4000 96 x 0.1ml
PSC - COMPATIBLE WITH: GE Healthcare MegaBACE 500 96 x 0.1ml
PSC - COMPATIBLE WITH: Hain-LIFESCIENCE GTQ-Cyclers 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Kyrathec SuperCyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: LongGene Scientific Instruments Co.
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: AgV
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-1 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-2 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-4 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAamplifier 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAquant 96 2plex 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAquant 96 5plex 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-4000 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-5000 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-Plus 96 x 0.2ml
PSC - COMPATIBLE WITH: Thermo Hybaid MultiBlock System 96 x 0.2ml
PSC - COMPATIBLE WITH: ThermoFisher Scientific Arktik Thermal Cyclers 96 x 0.2ml



PSC - COMPATIBLE WITH: VWR DuoCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR DuoCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR QuattroCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR QuattroCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR UnoCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR UnoCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TRobot 96 x 0.2ml
PSC - CAPACITY: 0.2ml
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%
PSC - GRIGLIA ALFANUMERICA: si