




copy of Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted


Cf 25pz - Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted. Cf 25pz



Product codes:

Reference: 12 893 0771

Product gallery:



Product description:

Piastra da 0.2ml semi skirted. Utilizzabile su sistemi PCR e PCR Real Time a profilo regolare. Questa piastra è inoltre perfettamente compatibile con sistemi ABI e Life Technologies

- Colore bianco
 - Regular Profile
 - Sub skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.1ml
PSC - FORM: semi-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V
PSC - CLOSING: through film



PSC - CLOSING: via caps
PSC - CLOSING: mat
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time FAST
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: transparent
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - PACKAGE: bulk
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: adjust
PSC - COMPATIBLE WITH: Agilent Technologies AriaDx Realtime PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Agilent Technologies AriaMx Realtime PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena Alpha SC 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler Monoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TOne 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TOne 96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Axygen® MaxyGene™ II 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer LifePro TC-96/G/H(b)A 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Opus Deepwell 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBLE WITH: Blue-Ray Biotech TurboCycler TCST-9610 96 x 0.2ml
PSC - COMPATIBLE WITH: Blue-Ray Biotech TurboCycler TCST-9620 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GenePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GeneTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio LifePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio LifeTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME GenePro Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME LifeExpress B-96D 96 x 0.2ml



PSC - COMPATIBLE WITH: eENZYME LifeExpress TC-96 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50a 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50i 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50l 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50s 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50t 384 Microplate
PSC - COMPATIBLE WITH: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ Maxi 96 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ MaxPro 96 x 0.2ml
PSC - COMPATIBLE WITH: Hain-LIFESCIENCE GTQ-Cycler 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Kyrtec SuperCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: LongGene Scientific Instruments Co.
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-1 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-2 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-4 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAamplifier 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAquant 96 2plex 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAquant 96 5plex 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR DuoCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR DuoCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR QuattroCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR QuattroCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR UnoCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR UnoCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems StepOnePlus™ Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX Duet Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml



PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories Opticon 2 96 x 0.2ml
PSC - COMPATIBLE WITH: Roche LightCycler ® 480 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Roche LightCycler ® 96 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO II 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO-thermoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: Acugen Systems AG-9600 thermal station 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena Alpha SC 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena Alpha SC 48-96 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler Twinblock 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler Twinblock Mix 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler2 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler2 48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler2 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 G 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 SG 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced Twin 48 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced Twin 48 G 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic Gradient XL96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic XL96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard Gradient SL 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard SL 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional TRIO 48 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional TRIO Combi 48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TRIO 48 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TRIO Combi 48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TRobot II 96 G 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TRobot II 96 SG 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC201 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC401 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC701 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC801 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 6 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 7 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems StepOne™ Real-Time PCR System 48 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems ViiA™ 7 Real-Time PCR System Fast 96 x 0.1ml
PSC - COMPATIBLE WITH: Barnstead Thermolyne Amplitron I 96 x 0.2ml
PSC - COMPATIBLE WITH: Barnstead Thermolyne Amplitron II 96 x 0.2ml



PSC - COMPATIBLE WITH: Bioer GeneExplorer GE-48D 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GeneExplorer GE-96G 96 x
PSC - COMPATIBLE WITH: 0.2ml
PSC - COMPATIBLE WITH: Bioer GenePro B-3048 30-48 30 x 0.5ml
PSC - COMPATIBLE WITH: Bioer GenePro B-3048 30-48 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GenePro B-48D 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer
PSC - COMPATIBLE WITH: GenePro B-96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GeneQ TC-24/H(b) 24 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GeneTouch™ Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GeneTouch™ Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Life ECO Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Life Express TC-48/H(t) 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Life Express TC-96/T/H(a) 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer LifeTouch TC-96/G/H(b)B 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Little Genius TC-24/H 24 x 0.2ml
PSC - COMPATIBLE WITH: Bioer XPcycler TC-XP-A 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer XPcycler TC-XP-D 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer XPcycler TC-XP-E 48-30 48 x 0.2ml
PSC - COMPATIBLE WITH: Bioer XPcycler TC-XP-G 96 x 0.2ml
PSC - COMPATIBLE WITH: Biomérieux GENE-UP 96 x 0.2ml Fast
PSC - COMPATIBLE WITH: Bioneer Exicycler 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Dual 48/48 Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch Dual 48/48 Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 6-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories MJ Mini 48-12 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC Tempo 96 Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 Dual 48/48 Fast
PSC - COMPATIBLE WITH: Biosan Ltd. TS-100 Thermo-Shaker 96 x 0.2ml
PSC - COMPATIBLE WITH: Boeco TC-PRO 3048 30-48 48 x 0.2ml
PSC - COMPATIBLE WITH: Boeco TC-PRO 48D 48 x 0.2ml
PSC - COMPATIBLE WITH: Boeco TC-PRO 96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GenePro™ Thermal Cyclers 30-48 48 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GenePro™ Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GeneQ™ Thermal Cyclers 24 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GeneTouch™ Thermal Cyclers 30-48 48 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GeneTouch™ Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBLE WITH: Chai Biotechnologies Inc. OpenQPCR Dual Channel 16 x 0.2ml
PSC - COMPATIBLE WITH: hai Biotechnologies Inc. OpenQPCR Single Channel 8 x 0.2ml
PSC - COMPATIBLE WITH: DuPont Qualicon BAX System Q7 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME GenePro Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME GeneQ Personal Thermal Cyclers 24 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME LifeExpress TC-48 48 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME MyGene MG-96 Gradient Thermal Cyclers 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME MyGene MG-96 Thermal Cyclers 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler gradient 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus eco 96 x 0.2ml



PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml
PSC - COMPATIBLE WITH: Ericom Deltacycler I 96 x 0.2ml
PSC - COMPATIBLE WITH: Ericom Deltacycler II 96 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ Maxi 48 x 0.2ml
PSC - COMPATIBLE WITH: sco Swift™ Maxi 48-30 48 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ Maxi 96 x 0.2ml gradient
PSC - COMPATIBLE WITH: Esco Swift™ MaxPro 48-30 48 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ MaxPro 48-48 48 x 0.2ml
PSC - COMPATIBLE WITH: Esco Swift™ MiniPro 24 x 0.2ml
PSC - COMPATIBLE WITH: EuroClone One Advanced 96 x 0.2ml
PSC - COMPATIBLE WITH: EuroClone One Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: EuroClone peqSTAR 96 HPL 96 x 0.2ml
PSC - COMPATIBLE WITH: EuroClone peqSTAR 96 ML 96 x 0.2ml
PSC - COMPATIBLE WITH: EuroClone peqSTAR 96 Universal 96 x 0.2ml
PSC - COMPATIBLE WITH: Finnzymes Instruments Piko™ 24-well 24-wells
PSC - COMPATIBLE WITH: Finnzymes Instruments TCA0096 Arktik 96 x 0.2ml
PSC - COMPATIBLE WITH: Finnzymes Instruments TCA4848 Arktik 2 x 48 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS1 48 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS1 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS1 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS2 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS4 48 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS4 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS4 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS4822 48 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GSX 48 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GSX 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GSX 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Hangzhou Jingle Scientific Instrument Co.
PSC - COMPATIBLE WITH: Helena Biosciences Phoenix Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Lab-Line Instruments Programmable Thermal Block 30-48 48 x 0.2ml
PSC - COMPATIBLE WITH: Lab-Line Instruments Programmable Thermal Block 96 x 0.2ml
PSC - COMPATIBLE WITH: Labnet MultiGene™ Mini Personal Thermal Cycler 24 x 0.2ml
PSC - COMPATIBLE WITH: Labnet MultiGene™ Thermal Cycler Gradient TC9600-G 96 x 0.2ml
PSC - COMPATIBLE WITH: Labnet MultiGene™ Thermal Cycler TC9600 96 x 0.2ml
PSC - COMPATIBLE WITH: Landgraf Varius M 96 x 0.2ml
PSC - COMPATIBLE WITH: Landgraf Varius V 80 16 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml AgV
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml



PSC - COMPATIBLE WITH: SensoQuest Labcycler Thermoblock 96
PSC - COMPATIBLE WITH: TaKaRa Bio Inc. PCR Thermal Cycler Dice TP950 96 x 0.1ml
PSC - COMPATIBLE WITH: TaKaRa Bio Inc. PCR Thermal Cycler Dice TP990 96 x 0.1ml
PSC - COMPATIBLE WITH: ThermoFisher Scientific Arktik Thermal Cycler 48 x 2ml
PSC - COMPATIBLE WITH: Whatman Biometra T3 Thermocycler 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T3 Thermocycler Combi 48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T3000 Thermocycler 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T3000 Thermocycler Combi 48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TOptical 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TPersonal 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TPersonal Combi 48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional TRIO 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional TRIO combi 48-30 48 x 0.2ml
PSC - COMPATIBLE WITH: Xian Tianlong Science & Technology DTC-4 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Xian Tianlong Science & Technology DTC-4G 96-77 96 x 0.2ml
PSC - CAPACITY: 0.1ml
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%
PSC - GRIGLIA ALFANUMERICA: si