




copy of Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted


Cf 25pz - Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted. Cf 25pz



Product codes:

Reference: 12 893 0801

Product gallery:



Product description:

Piastra da 0.2ml semi skirted. Utilizzabile su sistemi PCR e PCR Real Time a profilo regolare. Questa piastra è inoltre perfettamente compatibile con sistemi ABI e Life Technologies

- Colore bianco
 - Regular Profile
 - Sub skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.1ml
PSC - FORM: semi-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V



PSC - CLOSING: through film
PSC - CLOSING: via caps
PSC - CLOSING: mat
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time FAST
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: transparent
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - PACKAGE: bulk
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: adjust
PSC - COMPATIBLE WITH: Analytik-Jena Alpha SC 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler Monoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TOne 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TOne 96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Axygen® MaxyGene™ II 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GenePro B-96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer LifePro TC-96/G/H(b)A 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Line-Gene 9600 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Opus Deepwell 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBLE WITH: Blue-Ray Biotech TurboCycler TCST-9610 96 x 0.2ml
PSC - COMPATIBLE WITH: Blue-Ray Biotech TurboCycler TCST-9620 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GenePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio GeneTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio LifePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bulldog Bio LifeTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME GenePro Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME LifeExpress B-96D 96 x 0.2ml



PSC - COMPATIBLE WITH: eENZYME LifeExpress TC-96 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50a 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50i 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50l 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50s 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler X50t 384 Microplate
PSC - COMPATIBLE WITH: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBLE WITH: Hain-LIFESCIENCE GTQ-Cycler 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Kyrtec SuperCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: LongGene Scientific Instruments Co.
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-1 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-2 96 x 0.2ml
PSC - COMPATIBLE WITH: PCRMax AlphaCycler AC-4 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAamplifier 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAquant 96 2plex 96 x 0.2ml
PSC - COMPATIBLE WITH: Qiagen QIAquant 96 5plex 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-4000 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-5000 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-Plus 96 x 0.2ml
PSC - COMPATIBLE WITH: ThermoFisher Scientific Arktik Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR DuoCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR DuoCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR QuattroCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR QuattroCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR UnoCycler 96 x 0.2ml
PSC - COMPATIBLE WITH: VWR UnoCycler fast 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems StepOnePlus™ Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX Duet Real-Time PCR System 96 x 0.2ml



PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories Opticon 2 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. Opticon 2 System 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. Opticon System 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Roche LightCycler® 480 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Roche LightCycler® 96 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO II 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO-thermoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena FlexCycler2 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 G 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 SG 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic Gradient XL96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Basic XL96 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard Gradient SL 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TProfessional Standard SL 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC201 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC401 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC701 96 x 0.2ml
PSC - COMPATIBLE WITH: Apollo ATC801 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 6 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 7 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems ViiA™ 7 Real-Time PCR System Fast 96 x 0.1ml
PSC - COMPATIBLE WITH: Barnstead Thermolyne Amplitron I 96 x 0.2ml
PSC - COMPATIBLE WITH: Barnstead Thermolyne Amplitron II 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GeneTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer Life Express TC-96/T/H(a) 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer LifeTouch TC-96/G/H(b)B 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer XPcycler TC-XP-A 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC Tempo 96 Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: eENZYME MyGene MG-96 Gradient Thermal Cycler 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler gradient 96 x 0.2ml



PSC - COMPATIBLE WITH: Eppendorf Mastercycler gradient 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus eco 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml
PSC - COMPATIBLE WITH: Finnzymes Instruments TCA0096 Arktik 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS1 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS1 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS2 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS4 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GS4 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GSX 96 x 0.2ml
PSC - COMPATIBLE WITH: Gene Technologies Ltd. G-Storm GSX 96-48 96 x 0.2ml
PSC - COMPATIBLE WITH: SensoQuest Labcycler Thermoblock 96
PSC - COMPATIBLE WITH: TaKaRa Bio Inc. PCR Thermal Cyclers TP950 96 x 0.1ml
PSC - COMPATIBLE WITH: TaKaRa Bio Inc. PCR Thermal Cyclers TP990 96 x 0.1ml
PSC - COMPATIBLE WITH: Analytik-Jena qTower 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena qTOWER³ touch 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena SpeedCycler² 96 x 0.2ml
PSC - COMPATIBLE WITH: Analytik-Jena - Biometra TAdvanced 96 S 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Standard
PSC - COMPATIBLE WITH: Azure biosystems Cielo 3 Real Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Azure biosystems Cielo 6 Real Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bioer GeneExplorer GE-96G 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC Tempo Deepwell Thermal Cyclers 96-Deep Well Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: Biosystems MJ96+ 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler gradient 96-77 77 x 0.5ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus GSX1 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus GSX1e 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus SX1 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus SX1e 96 x 0.2ml
PSC - COMPATIBLE WITH: Landgraf Varius V 45 9 x 0.5ml
PSC - COMPATIBLE WITH: AgV MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Peqlab Primus 96 Advanced 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Roche cobas z 480 96 x 0.2ml
PSC - COMPATIBLE WITH: Seegene Inc. SEEAMP™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne Flexigene 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne Gene-E 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne Genius 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne Prime thermal cyclers 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-412 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne TC-512 96 x 0.2ml
PSC - COMPATIBLE WITH: Techne Touchgene Gradient 96 x 0.2ml
PSC - CAPACITY: 0.1ml
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%



Shop name

PSC - BINDING: proteins less 0.8%
PSC - GRIGLIA ALFANUMERICA: si