



Micropiastre PCR 48x0,2 Non-Skirted


Micropiastre PCR 48x0,2 neutre Non-skirted 50pz.



Product codes:

Reference: 12 893 0741

Product gallery:



Product description:

Micropiastre 48 x 0,2 ml Non-skirted
Parete sottile, profilo regolare, colore Neutro

Facilmente tagliabile in pezzi 16 o da 24 pozzetti.

50Pz

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.2ml
PSC - FORM: no-skirted
PSC - OPENING: normal
PSC - WELLS: 48
PSC - BOTTOM: a V
PSC - CLOSING: via caps
PSC - GRID: alphanumeric



PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time
PSC - APPLICATIONS: PCR Real-Time FAST
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: neutral
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - STACKABLE: you
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - MAIN FEATURES: PCR 48 wells for PCR and Real Time PCR no skirted cuttable
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: adjust
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml



PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TOptical-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
PSC - COMPATIBLE WITH: Stratagene-Robocycler Gradient 96-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-310-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-3100-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-3130-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-3700-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-3730-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-SimpliAmp™-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-Solid3-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml



PSC - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
PSC - CAPACITY: 0.2ml
PSC - PACKAGE: bulk
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: none
PSC - SIGNAL-STRENGTH: 100
PSC - INTERFERENCE SIGNAL: 100
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%