





Micropiastre PCR 96x0,1 neutre Non-skirted

Micropiastre PCR 96x0,1 neutre Non-skirted

Product codes:

Reference: 12 893 9061

Product gallery:



Product description:

Micropiastre in polipropilene Frosted PCR 96x0,1 no skirted, superficie qPCR
Low binding.

Disegnata per apparecchi BioRad tipo fast.
Tagliabile in piastre da 16/24/32 e 48 pozzetti.

Pareti uniformi e sottili.
Evaporazione minore del 3%

DNA, DNase, RNA, RNase e Pyrogen Free

Cf. 25pz



Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.1ml
PSC - FORM: no-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V
PSC - CLOSING: through film
PSC - CLOSING: via caps
PSC - CLOSING: mat
PSC - GRID: alphanumeric
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time
PSC - APPLICATIONS: PCR Real-Time FAST
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: neutral
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - STACKABLE: you
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - MAIN FEATURES: PCR for Real-Time and Real-Time FAST no skirted cuttable
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - OPTIONAL ACCESSORIES: the grid assembly white
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: PCR Real-Time FAST
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: low type FAST/0.1
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml



PSC - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-Opticon System-96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
PSC - COMPATIBLE WITH: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
PSC - CAPACITY: 0.2ml
PSC - PACKAGE: bulk
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%