



Strips 8 pozzetti qPCR da 0,1ml

Cf. 120 pezzi - Microtubi da 0,1ml per qPCR, profilo basso

Product codes:

Reference: 12 954 8391



Product gallery:



Product description:

Strips da 0,1ml, profilo basso, eccellenti per PCR Real-Time.

La scelta migliore per gli apparecchi ABI e BioRad per PCR Real-Time a basso profilo.

Tutta la serie di strips possiede una forma asimmetrica come riferimento orientativo, un'area scrivibile e dei fori per fissarli sulla griglia atta ad ospitare le strips per creare la propria piastra.

Cf. 120pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Leggenda:

Grado qPCR 1: sviluppato per segnale ottimo in PCR Real-Time

Grado qPCR 2: sviluppato per avere il segnale al top, esente da interferenze, in PCR Real-Time



Product features:

SP - TYPE: strip da 8 microtubi
SP - CAP: senza
SP - VOLUME: 0.1ml
SP - ABILITY: 0.2ml
SP - WELLS: 8
SP - PROFILE: low type FAST/0.1
SP - HEIGHT: 15.6 / 16.6mm ca
SP - SHAPE: strip sturdy
SP - FUND: tapered
SP - MATERIAL: PP
SP - APPLICATIONS: PCR
SP - APPLICATIONS: PCR Real Time
SP - APPLICATIONS: PCR Real Time FAST
SP - AUTOCLAVABLE: you
SP - COLOR: bianco
SP - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITY: no
SP - CONFORMITY: not subject to CE
SP - MAX TEMPERATURE: 100°
SP - MINIMUM TEMPERATURE: -30°C
THE SP - PACKAGE: envelope
SP - FEATURES: strip da 8x0.1 per PCR Real Time FAST e PCR extra robuste bianco
SP - OPTIONAL ACCESSORIES: the strip of caps
SP - OPTIONAL ACCESSORIES: griglia assemblaggio bianca
SP - SUITABLE FOR: PCR
SP - SUITABLE FOR: PCR Real Time
SP - SUITABLE FOR: PCR Real Time FAST
SP - SUITABLE FOR: storage
SP - WALLS: thin 0.30 mm
SP - TREATMENT SIGNAL: pareti bianche
SP - SIGNAL-STRENGTH: 180
SP - POSSIBLE INTERFERENCE SIGNAL: inferiore a 2
SP - RESISTANT: freezing at -20°C
SP - EVAPORATION SAMPLE: the lesser of 3%
SP - BINDING: DNA minor 0.3%
SP - BINDING: proteins less 0.8%
SP - COMPATIBLE WITH: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml



SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-Opticon System-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml



SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Roche-LightCycler® 480 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml