




tips gel, 200ul EC tips non-sterile in a sealed envelope

tips for gel EC-tips non-sterile - bags of 1,000 pieces-Cf.10000pz

**Product codes:**

Reference: 64 885 0251

Product gallery:**Product description:**

Ferrules made of transparent, with a section of the capillary is long and thin for gel loading allow extreme accuracy in sample loading in applications such as DNA sequencing and lelettroforesi of proteins. The capillary thin, it is particularly suitable for handling the gel in the wells is extremely narrow. Tips for gel available in bags of 1,000 pieces non-sterile and in rack with lid of 96 pieces sterile and free from DNA, Dnase and Rnase Apirogeni. Compatible with pipettes Gilson Pipetman® P20-P100-P200 Eppendorf® Biohit® Proline & mLine Socorex® Acura HTL® Discovery Nichiryo® Thermo Scientific Finnpiquette and the like. Also available with filter (barrier spray) to avoid contamination from the air and the cross-contamination of the samples. Material: polypropylene. Autoclavable.

Product features:

PNT - SIZE: nd
THE PNT - TYPE: gel loading
THE PNT - TYPE: without filter
PNT - CAPACITY MAX: 200ul
PNT - RANGE: 0.1-200ul
PNT - FILTER: absent
PNT - FORM: capillary
PNT - TIP: capillary round



PNT - STERILITY: no
PNT - PURITY: DNA / DNase / RNase / Pyrogen FREE
PNT - AUTOCLAVABLE: si
PNT - SIZED: no
PNT - COMPLIANCE: CE
PNT - APPLICATIONS: loading gel
PNT - MATERIAL: PP
PNT - COLOR: neutral
PNT - FEATURES: tip for gel loading by 200ul non-sterile in a sealed envelope
PNT - COMPATIBLE WITH: Biohit mLine
PNT - COMPATIBLE WITH: Biohit Proline
PNT - COMPATIBLE WITH: Eppendorf
PNT - COMPATIBLE WITH: Finnpipette
PNT - COMPATIBLE WITH: Socorex Acura
PNT - COMPATIBLE WITH: Nichiryo
PNT - COMPATIBLE WITH: Gilson Pipetman
PNT - COMPATIBLE WITH: HTL Disocvery
PNT - PACKAGING: in the envelope