




# tips for gel neutral from 0.1 to 200ul rack

tips for gel neutral from 0.1 to 200ul DNA Dnase Rnase free sterile apirogeni comp. GILSON EPPENDORF SOCOREX HTL NICHIRIO. rack 96pz-Cf.960pz

**Product codes:**

Reference: 64 885 9011

**Product gallery:****Product description:**

Ferrules made of transparent, with a section of the capillary is long and thin for gel loading allow extreme accuracy in sample loading in applications such as DNA sequencing and lelettroforesi of proteins. The capillary thin, it is particularly suitable for handling the gel in the wells is extremely narrow. Tips for gel available in bags of 1,000 pieces non-sterile and in rack with lid of 96 pieces sterile and free from DNA, Dnase and Rnase Apirogeni. Compatible with pipettes Gilson Pipetman® P20-P100-P200 Eppendorf® Biohit® Proline & mLine Socorex® Acura HTL® Discovery Nichiryo® Thermo Scientific Finnpiptette and the like. Also available with filter (barrier spray) to avoid contamination from the air and the cross-contamination of the samples. Material: polypropylene. Autoclavable.

**Product features:**

PNT - SIZE: nd  
THE PNT - TYPE: gel loading  
THE PNT - TYPE: without filter  
PNT - CAPACITY MAX: 200ul  
PNT - RANGE: 0.1-200ul  
PNT - FILTER: absent  
PNT - FORM: capillary  
PNT - TIP: capillary round  
PNT - STERILITY: no



---

PNT - PURITY: DNA / DNase / RNase / Pyrogen FREE  
PNT - AUTOCLAVABLE: si  
PNT - SIZED: no  
PNT - COMPLIANCE: CE  
PNT - APPLICATIONS: loading gel  
PNT - MATERIAL: PP  
PNT - COLOR: neutral  
PNT - FEATURES: tip for gel loading by 200ul sterile rack  
PNT - COMPATIBLE WITH: Biohit mLine  
PNT - COMPATIBLE WITH: Biohit Proline  
PNT - COMPATIBLE WITH: Eppendorf  
PNT - COMPATIBLE WITH: Finnpiquette  
PNT - COMPATIBLE WITH: Socorex Acura  
PNT - COMPATIBLE WITH: Nichiryo  
PNT - COMPATIBLE WITH: Gilson Pipetman  
PNT - COMPATIBLE WITH: HTL Disocvery  
PNT - PACKAGING: on the rack