



Microtubi qPCR 0.2ml con tappo piatto ottico - ultra trasp

Cf 1000pz - Microtubi qPCR 0.2ml con tappo piatto ottico - ultra trasp. Cf 1000pz

**Product codes:**

Reference: 12 879 6741

Product gallery:**Product description:**

Microtubi da 0.2ml con tappo ottico piatto anti-impronta. Compatibili con molti termociclatori per PCR e PCR Real Time che accettano plastica per biologia molecolare a profilo regolare. Come ABI / LT, Bio-Rad, Eppendorf e altri.

- Profilo regolare
- Ultra trasparenti
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- Tappo ottico
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

**Product features:**

MTB - CAP: dish
MTB - TYPE: PCR
MTB - VOLUME: 0.2ml
MTB - CAPACITY: 0.2ml
MTB - CLOSING: pressure
MTB - CAP FORABILE: nd
MTB - INSERTED CAP: no
MTB - OPTICAL AREA: 12.6mm
MTB - FORM: adjust
MTB - FUND: tapered
MTB - SELF STANDING: no
MTB - MATERIAL: PP Low Binding
MTB - SIZED: no
MTB - OPACITY: ultra-transparent
MTB - AREA WRITABLE: no
MTB - LABEL: absent
MTB - TEMPERATURE RANGE: -30°C a +100°C
MTB - TEMPERATURE MAX: 100°C
MTB - MINIMUM TEMPERATURE: -30°C
MTB - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - AUTOCLAVABLE: you
MTB - STERILITY: no
MTB - SIZE: h 20.9/21.8mm
MTB - FEATURES: microtubi da 0.2ml regular profile con tappo piatto per PCR Real-Time ultra-trasparenti - raccomandate quando il segnale è ottenuto dal fondo del microtubo
MTB - SUITABLE FOR: sample preparation
MTB - SUITABLE FOR: storage
MTB - SUITABLE FOR: PCR
MTB - SUITABLE FOR: PCR Real Time
MTB - SUITABLE FOR: molecular biology
MTB - RESISTANT: freezing at -20°C
MTB - PACKAGING: envelope bulk
MTB - WALLS: thin 0.3 mm
MTB - EVAPORATION SAMPLE: the lesser of 3%
MTB - BINDING: DNA minor 0.3%
MTB - BINDING: proteins less 0.8%
MTB - SIGNAL-STRENGTH: 100
MTB - INTERFERENCE SIGNAL: 100
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml



MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml



MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-SimpliAmp™-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
MTB - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml



MTB - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBLE WITH: Stratagene-Mx3000P™-96 x 0.2ml
MTB - COMPATIBLE WITH: Stratagene-Mx3005P™-96 x 0.2ml
MTB - COMPATIBLE WITH: Stratagene-Mx4000™-96 x 0.2ml
MTB - COMPATIBLE WITH: Stratagene-Robocycler Gradient 96-96 x 0.2ml
MTB - COMPATIBLE WITH: Whatman Biometra-TOptical-96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories iQ Real-Time PCR 96 x 0.2ml
MTB - COMPATIBLE WITH: Bio-Rad Laboratories iQ5 Real-Time PCR 96 x 0.2ml