



Piastra semi-skirted 96x0.1ml (LP) tagliabile - frosted

Cf 25pz - Piastra semi-skirted 96x0.1ml (LP) tagliabile - frosted. Cf 25pz

Product codes:

Reference: 12 893 8841



Product gallery:



Product description:

Piastra low profile da 0.1ml semi skirted. Utilizzabile su sistemi per PCR Real Time e PCR Real Time di tipo Fast.

La piastra può essere facilmente tagliata in sezioni da 3 file ed è spezzabile alle basse temperature (4°C)

- Colore neutro frosted
 - Low Profile
 - Semi skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.1ml
PSC - FORM: semi-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V



PSC - CLOSING: through film
PSC - CLOSING: via caps
PSC - CLOSING: mat
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time
PSC - APPLICATIONS: PCR Real-Time FAST
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: white
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - PACKAGE: bulk
PSC - MAIN FEATURES: piastra PCR 96 pozzetti da 0.1 per PCR e PCR Real Time e FAST semi-skirted
spezzabile a 4°C in tre strisce - FROSTED
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: PCR Real-Time FAST
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: low type FAST/0.1
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories Opticon 2 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories T100 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml



PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBLE WITH: Eppendorf Thermomixer comfort 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. Opticon 2 System 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. Opticon System 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml AgV
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO II 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO-thermoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch Dual 48/48 Fast
PSC - COMPATIBLE WITH: Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems StepOnePlus Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems StepOne Real-Time PCR System 48 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti 384-Well Thermal Cycler 384 microplate
PSC - COMPATIBLE WITH: Applied Biosystems Veriti 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti Dx 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Dual 48/48 Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories MiniOpticon 48-12 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories MJ Mini 48-12 48 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 Dual 48/48 Fast
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml



PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional TRIO combi 48-30 48 x 0.2ml
PSC - CAPACITY: 0.1ml
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%
PSC - GRIGLIA ALFANUMERICA: si