




Piastra semi-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted

Cf 25pz - Piastra semi-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted. Cf 25pz

Product codes:

Reference: 12 893 0441



Product gallery:



Product description:

Piastra da 0.2ml semi skirted. Utilizzabile su sistemi PCR e PCR Real Time a profilo regolare. Questa piastra è inoltre perfettamente compatibile con sistemi ABI e Life Technologies

- Colore neutro frosted
 - Regular Profile
 - Semi skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.2ml
PSC - FORM: semi-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V
PSC - CLOSING: through film
PSC - CLOSING: via caps



PSC - CLOSING: mat
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: PCR Real-Time
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: frosted
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - PACKAGE: bulk
PSC - MAIN FEATURES: piastra PCR 96 pozzetti da 0.2 per PCR e PCR Real Time semi-skirted - FROSTED
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: adjust
PSC - COMPATIBLE WITH: Applied Biosystems 310 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3100 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3130 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3730 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Standard
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 2700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 2720 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 Dual 96-Well 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 Aluminum 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 Gold 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems ProFlex™ PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems SimpliAmp™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems Solid3 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ 96-Well Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ Dx 96-Well Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems ViiA™ 7 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml



PSC - COMPATIBLE WITH: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories iCycler 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories iQ Real-Time PCR 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories iQ5 Real-Time PCR 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories MyCycler 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories MyiQ Real-Time PCR 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories Opticon 2 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories PTC-100 96 x 0.2ml

PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Deep Well

PSC - COMPATIBLE WITH: Bio-Rad Laboratories S1000 96-Well Fast

PSC - COMPATIBLE WITH: Bio-Rad Laboratories T100 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler 96-77 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler gradient 96-77 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus eco 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus gradient 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus GSX1 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus GSX1e 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus SX1 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler nexus SX1e 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Mastercycler pro S 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Thermomixer C 96 x 0.2ml

PSC - COMPATIBLE WITH: Eppendorf Thermomixer comfort 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-100 96 x 0.2ml AgV

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml

PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml



PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBLE WITH: Stratagene Mx3000P™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Stratagene Mx3005P™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TOptical 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO II 96 x 0.2ml
PSC - COMPATIBLE WITH: Whatman Biometra UNO-thermoblock 96 x 0.2ml
PSC - COMPATIBLE WITH: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7300 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7500 Real-Time PCR System 96 x 0.2ml Standard
PSC - COMPATIBLE WITH: Applied Biosystems ABI PRISM 7000 SDS 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems ABI PRISM 7700 SDS 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 5700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems MiniAmp Plus Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems MiniAmp Thermal Cycler 96 x 0.2ml
PSC - CAPACITY: 0.2ml
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%
PSC - GRIGLIA ALFANUMERICA: si