




Strip 8x0.1ml (LP) con tappo ottico attaccato - pareti frosted

Cf 120pz - Strip 8x0.1ml (LP) con tappo ottico attaccato - pareti frosted. Cf 120pz



Product codes:

Reference: 12 954 8091

Product gallery:



Product description:

Strip da 8 microtubi con tappo ottico attaccato compatibile con molti strumenti che accettano plastica per biologia molecolare di tipo FAST - Low Profile

- Disegnata per applicazioni in PCR e PCR Real Time è ideale inoltre per sistemi come ABI / Life technologies di tipo Fast, Roche, Bio Rad e altri. Il tappo attaccato ad ogni pozzetto possiede un'ampia area ottica di oltre 12mm

- Tappo attaccato su ogni singolo pozzetto

- Versione Low Profile da 0.1ml

- Colore neutro frosted

- Pareti sottili

- Miscela in PP specifico per applicazioni in PCR

- Basso grado di evaporazione

- Bassa adesione di DNA e proteine sulle pareti del tubo

- Tappo ottico

- DNA, RNA, DNase, RNase e Pyrogen FREE

- Autoclavabili

Product features:

SP - TYPE: strip da 8 microtubi



SP - CAP: attaccato
SP - CAP: piano
SP - CAP: ottico
SP - VOLUME: 0.2ml
SP - ABILITY: 0.2ml
SP - WELLS: 8
SP - PROFILE: low type FAST/0.1
SP - HEIGHT: 15.6 / 16.6mm ca
SP - SHAPE: striscia sottile
SP - FUND: tapered
SP - MATERIAL: PP
SP - APPLICATIONS: PCR
SP - APPLICATIONS: PCR Real Time
SP - APPLICATIONS: PCR Real Time FAST
SP - AUTOCLAVABLE: you
SP - COLOR: neutral
SP - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITY: no
SP - CONFORMITY: not subject to CE
SP - MAX TEMPERATURE: 100°C
SP - MINIMUM TEMPERATURE: -30°C
THE SP - PACKAGE: envelope
SP - FEATURES: strip da 8x0.1 per PCR Real Time FAST e PCR a striscia sottile con tappo attaccato - frosted
SP - OPTIONAL ACCESSORIES: griglia assemblaggio bianca
SP - OPTIONAL ACCESSORIES: griglia assemblaggio gialla
SP - SUITABLE FOR: PCR
SP - SUITABLE FOR: PCR Real Time
SP - SUITABLE FOR: PCR Real Time FAST
SP - WALLS: thin 0.30 mm
SP - TREATMENT SIGNAL: walls of frosted
SP - SIGNAL-STRENGTH: 125
SP - POSSIBLE INTERFERENCE SIGNAL: 70
SP - RESISTANT: freezing at -20°C
SP - EVAPORATION SAMPLE: the lesser of 3%
SP - BINDING: DNA minor 0.3%
SP - BINDING: proteins less 0.8%
SP - COMPATIBLE WITH: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
SP - COMPATIBLE WITH: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep



Well

SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-Opticon System-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE WITH: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml



SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Roche-LightCycler® 480 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml