



Strip 8x0.2ml (RP) extra-rigide - Ultra-trasp

Cf 120pz - Strip 8x0.2ml (RP) extra-rigide - Ultra-trasp. Cf 120pz

Product codes:

Reference: 12 954 3011

Product gallery:



Product description:

Strip da 8 microtubi compatibile con molti strumenti per PCR e PCR Real Time che accettano plastica per biologia molecolare a profilo regolare e alcuni modelli di sequenziatori

- Disegnata per applicazioni in PCR e PCR Real Time , compatibile anche con apparecchi come ABI / Life Technologies, Bio-Rad, Eppendorf e altri
- Striscia robusta non flessibile
- Versione profilo regolare da 0.2ml
- Colore neutro ultra trasparente
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

Product features:

- SP - TYPE: strip da 8 microtubi
- SP - CAP: senza
- SP - VOLUME: 0.2ml
- SP - ABILITY: 0.2ml
- SP - WELLS: 8
- SP - PROFILE: adjust



SP - HEIGHT: 20.9 / 21.8mm ca
SP - SHAPE: strip sturdy
SP - FUND: tapered
SP - MATERIAL: PP
SP - APPLICATIONS: PCR
SP - APPLICATIONS: PCR Real Time
SP - AUTOCLAVABLE: you
SP - COLOR: neutral
SP - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITY: no
SP - CONFORMITY: not subject to CE
SP - MAX TEMPERATURE: 100°C
SP - MINIMUM TEMPERATURE: -30°C
THE SP - PACKAGE: envelope
SP - FEATURES: strip da 8x0.2 per PCR Real Time e PCR a striscia robusta - ultra trasparenti
SP - OPTIONAL ACCESSORIES: griglia assemblaggio gialla
SP - SUITABLE FOR: PCR
SP - SUITABLE FOR: PCR Real Time
SP - WALLS: thin 0.30 mm
SP - TREATMENT SIGNAL: nessuno
SP - SIGNAL-STRENGTH: 100
SP - POSSIBLE INTERFERENCE SIGNAL: 100
SP - RESISTANT: freezing at -20°C
SP - EVAPORATION SAMPLE: the lesser of 3%
SP - BINDING: DNA minor 0.3%
SP - BINDING: proteins less 0.8%
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml



SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-310-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3100-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3130-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3500 Genetic Analyzer-8-capillary
SP - COMPATIBLE WITH: Applied Biosystems-3700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3730-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-SimpliAmp™-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-Solid3-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml



SP - COMPATIBLE WITH: Stratagene-Robocycler Gradient 96-96 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TOptical-96 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE WITH: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
SP - COMPATIBLE WITH: Bio-Rad Laboratories-iCycler-48 x 0.2ml
SP - COMPATIBLE WITH: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-3500xL Genetic Analyzer-24-capillary
SP - COMPATIBLE WITH: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-MiniAmp Plus Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE WITH: Applied Biosystems-MiniAmp Thermal Cycler-96 x 0.2ml