



Agitatore orbitale Orbi-Shaker XL da 30 a 300 rpm orbita 19mm. Capacità 30kg

Agitatore orbitale Orbi-Shaker XL da 30 a 300 rpm orbita 19mm. Capacità carico 30kg - Completo di piattaforma con tappetino antiscivolo 450x450mm



Product codes:

Reference: 94 907 0012

Product gallery:



Product description:

Agitatore orbitale Orbi-Shaker XL da 30 a 300 rpm orbita 19mm.

Modello da banco ad alta capacità

Orbita di 19mm per ottimale miscelazione ed areazione dei campioni

Cambio istantaneo delle clamp per beute grazie alla piattaforma Magic Clamp e relativi accessori -acquistabili separatamente -

Capacità carico 30kg

- Completo di piattaforma con tappetino antiscivolo 450x450mm

Product features:

STG - TYPE: orbital shaker

STG - CAPACITY: contenitori vario tipo

STG - CAPACITY: 6 micropiastre

STG - FLOW: 30kg

STG - PLATFORM DIMENSIONS CM.: 45x45

STG - MOVEMENT: orbital

STG - ORBIT: 19mm

STG - RPM MAXIMUM: 300

STG - RPM MINIMUM: 30

STG - INCREASES RPM: 1

STG - MODE OPERATIONS: timed / manual

STG - INTERFACE: display



STG - DISPLAY: LED

STG - OPERATING CONDITIONS AMB.: da +4°C a 45°C

STG - FEATURES: agitatore alta capacità con tappetino magnetico MagicClamp per istantaneo montaggio accessori

TSG - SIZE: 47x54x12.5cm

STG - WEIGHT: 29kg

STG - POWER: 230V - 50-60Hz

STG - ACCESSORIES INCLUDED: standard platform

STG - OPTIONAL ACCESSORIES: platforms

STG - OPTIONAL ACCESSORIES: piattaforma MagicClamp

STG - CAN BE USED WITH: the flasks 500ml

STG - CAN BE USED WITH: bottiglie 500ml

STG - CAN BE USED WITH: provette Ø14mm

STG - CAN BE USED WITH: beute 25ml

STG - CAN BE USED WITH: beute 50ml

STG - CAN BE USED WITH: beute 125ml

STG - CAN BE USED WITH: beute 3000ml

STG - CAN BE USED WITH: buste / sacche

STG - CAN BE USED WITH: tubes

STG - CAN BE USED WITH: the flasks 2000ml

STG - CAN BE USED WITH: the flasks 250ml

STG - CAN BE USED WITH: the flasks 1000ml

STG - CAN BE USED WITH: vessels low

STG - CAN BE USED WITH: plates multiwell

STG - CAN BE USED WITH: tubes 50 ml type Falcon

STG - CAN BE USED WITH: tubes 15ml type Falcon

STG - CAN BE USED WITH: tubes Ø30mm

STG - CAN BE USED WITH: microtubes 2.0 ml

STG - CAN BE USED WITH: microtubes 1.5 ml