



# Termoblocco MyBlock Mini da 20° a 100°C a blocchi intercambiabili da 0.2 a 50ml

Termoblocco MyBlock Mini da +20° (+5°C amb) a 100°C a blocchi intercambiabili da 0.2 a 50ml. Con coperchio trasparente. Senza blocchi da acquistare separatamente

## Product codes:

Reference: 94 912 0120

## Product gallery:



## Product description:

Con la sua struttura compatta, l'ampia gamma di temperature e le opzioni di più blocchi intercambiabili, myBlock Mini è lo strumento perfetto per le incubazioni di laboratorio. I semplici controlli tramite touchpad e un display digitale consentono la selezione della temperatura "set and walk away" con un'accuratezza eccellente.

Perfetto per l'uso da banco, l'adattatore opzionale da 12 V espande l'utilità di myBlock Mini all'utilizzo. Con meno di 11cm di larghezza, myBlock Mini è davvero il primo incubatore a blocchi intercambiabili personale.

## Product features:

TBH - TYPE: heater  
TBH - No BLOCKS: 1  
TBH - WITH COOLING: no  
TBH - TEMPERATURE MAXIMUM: +100°C  
TBH - MINIMUM TEMPERATURE SET: 20°C  
TBH - THE MINIMUM TEMPERATURE REACHED: +5°C above ambient  
TBH - POWER HEATING: 50W  
TBH - RESOLUTION SETTING: 0,1°C  
TBH - PRECISION: 0,5°C



---

TBH - UNIFORMITY: 0,2°C  
TBH - DISPLAY: LED  
TBH - TIMER: -  
TBH - PROGRAMS: not provided  
TBH - AUTONOMY: -  
TBH - COMPLIANCE: CE  
TBH - FREE: no  
TBH - INDEPENDENT BLOCK: no  
TBH - OPTIONAL ACCESSORIES: adattatore 12V  
TBH - OPTIONAL ACCESSORIES: blocchi varie misure  
TBH - POWER: 240 V 50/60 Hz  
TBH - WARRANTY: 12 months on manufacturing defects  
TBH - SIZE: 102x150x115 mm  
TBH - WEIGHT: 1 kg  
TBH - FEATURES: termoblocco compatto con temperature a blocchi intercambiabili estremamente versatile  
TBH - TO BE USED WITH: the cuvette  
TBH - TO BE USED WITH: vials 12x32  
TBH - TO BE USED WITH: provette Ø30mm  
TBH - TO BE USED WITH: provette 50ml  
TBH - TO BE USED WITH: provette 15ml  
TBH - TO BE USED WITH: microtubi 5.0ml  
TBH - TO BE USED WITH: provette Ø13mm  
TBH - TO BE USED WITH: provette Ø12mm  
TBH - TO BE USED WITH: strip 8x0.2ml  
TBH - TO BE USED WITH: microtubes 2.0 ml  
TBH - TO BE USED WITH: microtubes 1.5 ml  
TBH - TO BE USED WITH: microtubes 0.5 ml  
TBH - TO BE USED WITH: microtubes 0.2 ml  
TBH - IDEAL FOR: incubazione  
TBH - UTILIZABILE CON: provette Ø16mm