



# Agarosio LM Low Melting Point. Da 200bp a 25kb. Cf 100g

Agarosio LM Low Melting Point per app in gel processing. Ideale per frammenti da 200bp a 25kb. Strength 200g/cm (1%) / M.T. 65°C / EEO 0.1. RNase / Dnase e protease FREE. Cf 100g



## Product codes:

Reference: 94 896 0213

## Product gallery:



## Product description:

Agarosio LM -Low Melting point -

Agarosio altamente purificato e scelta perfetta quando si vogliono effettuare manipolazioni di DNA in-gel.

Raffinato utilizzato un processo avanzato che esclude l'utilizzo di solventi organici, rendendo il prodotto finito puro e con impatto ambientale significativamente ridotto.

Questo agarosio è ideale per la separazione di frammenti da 200bp a 25kb.

Perfetto per la preparazione di gel con basso punto di fusione a temperature inferiori ai 65°C. L'agarosio rimane liquido a 37°C permettendo procedure come restrizione e labeling all'interno del gel.

Il prodotto permette la creazione di un gel trasparente, con bassa EEO e basso background. RNase, DNase, Protease and Endonuclease FREE

## Product features:

RGT - TYPE: agarosio

RGT - COMPLIANCE: LM



---

RGT - FORM: dust  
RGT - FRAGMENTS: 200 - 25kbp  
RGT - POINT FUSION: 65°C  
RGT - POINT GELATION: 26-30°C  
RGT - APPLICATION CLONING: 200  
RGT - APPLICATIONS WITH BACTERIA / PATHOGENS: 25.000  
RGT - APPLICATIONS WITH THE DNA: 200g/cm  
RGT - APPLICATIONS WITH RESTRICTION ENZYMES: 0.1  
RGT - APPLICATIONS WITH LIPIDS: 10%  
RGT - APPLICATIONS WITH RNA: no detected  
RGT - MARKING APPLICATIONS: no detected  
RGT - APPLICATIONS IN MOLECULAR BIOLOGY: no detected  
RGT - APPLICATIONS IN CHEMISTRY: barattolo  
RGT - APPLICATIONS PROTEINS / ENZYMES: processi in-gel  
RGT - APPLICATIONS PROTEINS / ENZYMES: recupero frammenti  
RGT - APPLICATIONS PROTEINS / ENZYMES: separazione frammenti DNA  
RGT - APPLICATION REVERSE TRANSCRIPTION: for research use only  
RGT - FEATURES: gel di agarosio Low Melting per applicazioni che prevedono il recupero dopo corsa elettroforetica o applicazioni in-gel  
RGT - STORAGE: environment  
RGT - TRANSPORT: ambient temperature  
RGT- PROTEASE: no detected