

## Additive Screen 02 (A2), 1.0 ml

0.1 M Cadmium chloride hydrate,

, Individual reagents from Additive Screen and Additive Screen HT for screening, optimization, and production.

**Product codes:** 

Reference: HR 242802



**Product gallery:** 



## **Product description:**

- 1 ml
- Formulated in Type 1+ ultrapure water: 18.2 megaohm-cm resistivity at 25°C, 5 ppb Total Organic Carbon, bacteria free (1 Bacteria (CFU/ml)), pyrogen free (0.03 Endotoxin (EU/ml)), RNase-free (0.01 ng/mL) and DNase-free (4 pg/uL)

HR2-428-\*\*

\*\* Refers to the reagent number in the kit. For example, reagent number 1 = HR2-428-01

The individual reagents 1 (A1) to 96 (H12) in the Additive Screen and Additive Screen HT are identical and use the same Additive Screen Custom Shop catalog number. Reagent is supplied in a 1.0 ml volume in an o-ring screw cap polypropylene tube.

June 2008 reagent 89(H5) changed to 40% v/v Formamide.

For Additive Screen reagent 51(E3) Dextran sulfate, please order HR2-428-51A for Additive Screen kits lot 2428xx where xx = 11-12, and order HR2-428-51B for Additive Screen kits lot 2428zz where zz = 1-10 or >12.



For Additive Screen HT reagent E3(51) Dextran sulfate, please order HR2-428-51A for Additive Screen HT kits lots 213806, 213807 and 213808, and order HR2-428-51B for Additive Screen HT kits lot 2138zz where zz = 1-5 or >8.

For Additive Screen reagent 52 (E4), please order HR2-428-52A for Additive Screen kits Lot 2428xx where xx ? 19 contain 40% v/v Pentaerythritol ethoxylate (3/4 EO/OH) and order HR2-428-52B for Additive Screen kits Lot 2428yy where yy ? 20 contain 40% v/v Polyethylene glycol 300 as reagent 52 (E4).

For Additive Screen HT reagent 52 (E4), please order HR2-428-52A for Additive Screen kits Lot 2138xx where xx ? 19 contain 40% v/v Pentaerythritol ethoxylate (3/4 EO/OH) and order HR2-428-52B for Additive Screen kits Lot 2138yy where yy ? 20 contain 40% v/v Polyethylene glycol 300 as reagent 52 (E4).

For Additive Screen reagent 68(F8) CYMAL-7, please order HR2-428-68A for Additive Screen kits lot 2428xx where xx = 07-17, and order HR2-428-68B 1,2,3-Heptanetriol for Additive Screen kits lot 2428zz where zz is greater than or equal to 18.

For Additive Screen HT reagent F8(68) CYMAL-7, please order HR2-428-68A for Additive Screen HT kits lots 2138xx where xx = 03 through 15, and order HR2-428-68B 1,2,3-Heptanetriol for Additive Screen HT kits lot 2138zz where zz is greater than or equal to 16.

Per maggiori informazioni visita il sito https://hamptonresearch.com/

## **Product features:**

CRF - TIPO: Individual Additive Screen Reagents