



Strips 8 tappi ottici piani anti-impronta

Cf. 120 strips da 8 tappi - Strips di 8 tappi ottici piani anti-impronta

Product codes:

Reference: 12 864 8451



Product gallery:



Product description:

Strips di 8 tappi ottici piani anti-impronta

La superficie rientrante evita che l'area ottica venga "contaminata" da impronte accidentali

Utilizzabili in tutte le procedure in PCR Real Time

Pareti sottili ed uniformi, ottimi per chiudere le piastre e le strips per le quali sono indicate come accessorio.

Ottimizzate per PCR Real Time

Cf. 120 strips da 8 tappi



Product features:

FLM - TIPO DE: la tira de tapas
FLM - COMPATIBILIDAD: strip PCR
FLM - COMPATIBILIDAD: placas de PCR
FLM - COMPATIBILIDAD: placas de PCR en Tiempo Real
FLM - COMPATIBILIDAD: strip PCR Real-Time
FLM - MATERIAL: PP
FLM - POZOS: 8
FLM - PUREZA: DNA / RNA / DNase / Rnase / Pyrogen FREE
FLM - ESTERILIDAD: no
FLM - FORMA: reentrada
FLM - FORMA: tira resistente tapas
FLM - TEMPERATURA MAX: 100°C
FLM - TEMPERATURA MÍNIMA: -30°C
FLM - ÁREA ÓPTICA: 12.6mm
FLM - ÁREA DE ESCRITURA: usted
FLM - ESTERILIZABLE EN AUTOCLAVE: usted
FLM - CODIFICACIÓN DE COLOR: no
FLM - ETIQUETA: no
FLM - CUMPLIMIENTO: no sujeto a CE
FLM - CARACTERÍSTICAS: strip da 8 tappi piani anti-interferenza extra-robusti
FLM - APROPIADO PARA: biología molecular
FLM - APROPIADO PARA: PCR
FLM - APROPIADO PARA: PCR Real Time
FLM - APROPIADO PARA: PCR Real Time FAST
FLM - APROPIADO PARA: almacenamiento
FLM - APROPIADO PARA: almacenamiento a -20°C
FLM - COMPATIBLE CON: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBLE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBLE CON: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml



FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
FLM - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-SimpliAmp™-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-Solid3-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
FLM - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
FLM - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
FLM - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
FLM - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
FLM - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
FLM - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
FLM - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml



FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
FLM - COMPATIBLE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-Opticon System-96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBLE CON: Roche-cobas z 480-96 x 0.2ml
FLM - COMPATIBLE CON: Roche-LightCycler ® 480 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE CON: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE CON: Stratagene-Mx3000P™-96 x 0.2ml
FLM - COMPATIBLE CON: Stratagene-Mx3005P™-96 x 0.2ml
FLM - COMPATIBLE CON: Stratagene-Mx4000™-96 x 0.2ml
FLM - COMPATIBLE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TOptical-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
FLM - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast



FLM - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: Roche-COBAS Taqman 48 Analyzer-24 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TPersonal-48 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
FLM - COMPATIBLE CON: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
FLM - COMPATIBLE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
FLM - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE CON: Roche-LightCycler® Nano-32 x 0.2ml
FLM - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml