



Microtubi PCR da 0,5ml tappo piatto, graduati

Cf. 1000 pezzi - Microtubi da 0,5ml per PCR, tappo piatto, graduati

Product codes:

Reference: 12 879 3051



Product gallery:



Product description:

Microtubi da 0,5ml, tappo piatto, graduati.

Compatibili con molti apparecchi per PCR che accettano provettine da 0,5ml con tappo piatto.

Profilo regolare.

Cf. 1000pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Product features:

MTB - CAP: plato

MTB - CUMPLIMIENTO: no sujeto a CE

MTB - TIPO DE: PCR



MTB - TIPO DE: Eppendorf
MTB - VOLUMEN: 0.5ml
MTB - CAPACIDAD: 0.5ml
MTB - CAP FORABILE: no
MTB - INSERTA LA TAPA: no
MTB - ÁREA ÓPTICA: 5mm
MTB - FORMULARIO: ajustar
MTB - FONDO: cónico
MTB - AUTO DE PIE: no
MTB - MATERIAL: PP
MTB - TAMAÑO: usted
MTB - OPACIDAD: ultra-transparente
MTB - ÁREA DE ESCRITURA: no
MTB - RANGO DE TEMPERATURA: -30°C a +100°C
MTB - TEMPERATURA MAX: 100°C
MTB - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - ESTERILIZABLE EN AUTOCLAVE: usted
MTB - ESTERILIDAD: no
MTB - TAMAÑO: h 30mm ca
MTB - CARACTERÍSTICAS: microtubi da 0.5ml con tappo piatto graduati per PCR e PCR Real-Time
MTB - APROPIADO PARA: preparación de la muestra
MTB - APROPIADO PARA: almacenamiento
MTB - APROPIADO PARA: PCR
MTB - APROPIADO PARA: PCR Real Time
MTB - RESISTENTE: la congelación a -20°C
MTB - EMBALAJE: sobres a granel
MTB - PAREDES: fina de 0,3 mm
MTB - UNIÓN: las proteínas de menos de 0,8%
MTB - UNIÓN: Menor del ADN 0.3%
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-48 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 77 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler-20 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler-20 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TGradient-48 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TPersonal-20 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TPersonal Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional-60 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional Standard-60 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO-30 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-18 18 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TRIO Thermoblock-20 x 0.5ml



MTB - COMPATIBLE CON: Whatman Biometra-UNO II-48 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-40 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-GeneCycler-12 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 480-48 x 0.5ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-60 x 0.5ml
MTB - COMPATIBLE CON: Applied Biosystems-Veriti™ 60-Well Thermal Cycler-60 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-60 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-MJ Mini-48-12 12 x 0.5ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-60 x 0.5ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 77 x 0.5ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96-77 77 x 0.5ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler Personal-25-16 16 x 0.5ml
MTB - COMPATIBLE CON: Eppendorf-Thermomixer C-24 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-100-60 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-150 MiniCycler-16 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-60 Alpha unit 60 x 0.5ml