



Piastra no skirted, prof. regolare, tagliabile. Tipo universale


Cf. 25 pezzi - Piastre no-skirted, profilo regolare, tagliabili, tipo universale

Product codes:

Reference: 12 893 4121



Product gallery:



Product description:

Piastre da 0,2ml, profilo regolare, eccellenti per PCR Real-Time. Colore neutro

Compatibilità con la maggior parte degli apparecchi per PCR e qPCR a profilo regolare. Tagliabili.

Cf. 25pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Leggenda:

Grado qPCR 1: sviluppato per segnale ottimo in PCR Real-Time

Grado qPCR 2: sviluppato per avere il segnale al top, esente da interferenze, in PCR Real-Time



Product features:

PSC - TIPO DE: placas de PCR
PSC - VOLUMEN: 0.2ml
PSC - FORMULARIO: no-skirted
PSC - APERTURA: normal
PSC - POZOS: 96
PSC - PARTE INFERIOR: a V
PSC - CIERRE: a través del cine
PSC - CIERRE: a través de las tapas de
PSC - CIERRE: mat
PSC - GRID: alfanumérico
PSC - APLICACIONES: PCR
PSC - APLICACIONES: PCR Real-Time
PSC - APLICACIONES: biología molecular
PSC - ESTERILIZABLE EN AUTOCLAVE: uested
PSC - COLOR: neutral
PSC - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - ESTERILIDAD: no
PSC - MATERIAL: PP
PSC - CUMPLIMIENTO: no sujeto a CE
PSC - APILABLE: uested
PSC - TEMPERATURA MAX: 100°C
PSC - TEMPERATURA MÍNIMA: 30°C
PSC - CARACTERÍSTICAS PRINCIPALES: placa de PCR de 96 pocillos de PCR y PCR en Tiempo Real no se puede cortar con un faldón
PSC - ACCESORIOS OPCIONALES: la tira de tapas
PSC - ACCESORIOS OPCIONALES: la película de los selladores
PSC - ACCESORIOS OPCIONALES: esteras de selladores
PSC - APROPIADO PARA: PCR
PSC - APROPIADO PARA: PCR Real-Time
PSC - APROPIADO PARA: almacenamiento
PSC - APROPIADO PARA: preparación de la muestra
PSC - PERFIL: ajustar
PSC - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
PSC - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
PSC - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
PSC - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
PSC - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
PSC - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
PSC - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
PSC - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
PSC - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml



PSC - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TOptical-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
PSC - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
PSC - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
PSC - COMPATIBLE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-310-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-3100-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-3130-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-3700-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-3730-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml



PSC - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-SimpliAmp™-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
PSC - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
PSC - CAPACIDAD: 0.2ml
PSC - PAQUETE: bulk
PSC - PAREDES: delgada de 0.35 mm
PSC - TRATAMIENTO DE LA SEÑAL: ninguno
PSC - INTENSIDAD DE SEÑAL: 100
PSC - INTERFERENCIA DE LA SEÑAL: 100
PSC - RESISTENTE: la congelación a -20°C
PSC - EVAPORACIÓN DE LA MUESTRA: el menor de 3%
PSC - UNIÓN: Menor del ADN 0.3%
PSC - UNIÓN: las proteínas de menos de 0,8%