



Strips in combipack 8x0,2 con tappi a cupola attaccati

Strips PCR EU 8x0,2 con tappi a cupola attaccati. Conf. da 120 strips.

Product codes:

Reference: 12 954 2551

Product gallery:



Product description:

Strips da 8x0,2 in Combipack con tappi a cupola attaccati.

In PP opacizzato, Colore Neutro

Progettato per le applicazioni PCR-Real Time

DNA, RNA, DN-RNase e Pyrogene Free.

In polipropilene low binding

Conf. Busta da 120pz

Product features:

SP - TIPO DE: strip da 8 microtubi

SP - CAP: attaccato



SP - CAP: semi-cupola
SP - VOLUMEN: 0.2ml
SP - CAPACIDAD: 0.2ml
SP - POZOS: 8
SP - PERFIL: ajustar
SP - ALTURA: 20.9 / 21.8mm ca
SP - FORMA: striscia sottile
SP - FONDO: cónico
SP - MATERIAL: PP
SP - APLICACIONES: PCR
SP - APLICACIONES: PCR Real Time
SP - ESTERILIZABLE EN AUTOCLAVE: uested
SP - COLOR: neutral
SP - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - ESTERILIDAD: no
SP - CONFORMIDAD: no sujeto a CE
SP - TEMPERATURA MAX: 100°
SP - TEMPERATURA MÍNIMA: -30°C
EL SP - PAQUETE: sobres
SP - CARACTERÍSTICAS: strip da 8x0.2 per PCR Real Time e PCR con tappo a semi-cupola attaccato ultra-trasparenti
SP - APROPIADO PARA: PCR
SP - APROPIADO PARA: PCR Real Time
SP - APROPIADO PARA: almacenamiento
SP - PAREDES: delgado 0,30 mm
SP - TRATAMIENTO DE LA SEÑAL: nessuno
SP - INTENSIDAD DE SEÑAL: 100
SP - A LA POSIBLE INTERFERENCIA DE LA SEÑAL: 100
SP - RESISTENTE: la congelación a -20°C
SP - EVAPORACIÓN DE LA MUESTRA: el menor de 3%
SP - UNIÓN: Menor del ADN 0.3%
SP - UNIÓN: las proteínas de menos de 0,8%
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml



SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml



SP - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml