



Strip di tubi Low-profile

Q.Tà DISPONIBILE: Ultimi Pezzi - Vedi Descrizione

Product codes:

Reference: 12 954 1191

Product gallery:



Product description:

Strips da 8/12 provette da 0,2ml Low-Profile, certificate per PCR.
Colore Naturale

Disponibilità:

Da 8 provette: 11 Conf. da 120 strips

Da 12 provette: 14 Conf. da 80 strips

Product features:

SP - TIPO DE: tira de 12 microtubos

SP - CAP: senza

SP - VOLUMEN: 0.1ml

SP - CAPACIDAD: 0.2ml

SP - POZOS: 12

SP - PERFIL: tipo de baja RÁPIDO/0.1

SP - ALTURA: 15.6 / 16.6mm ca



SP - FORMA: tira resistente
SP - FONDO: cónico
SP - MATERIAL: PP
SP - APLICACIONES: PCR
SP - APLICACIONES: PCR Real Time
SP - APLICACIONES: PCR Real Time FAST
SP - ESTERILIZABLE EN AUTOCLAVE: usted
SP - COLOR: neutral
SP - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - ESTERILIDAD: no
SP - CONFORMIDAD: no sujeto a CE
SP - TEMPERATURA MAX: 100°
SP - TEMPERATURA MÍNIMA: -30°C
EL SP - PAQUETE: sobres
SP - CARACTERÍSTICAS: tira de 12x0.1 por PCR en Tiempo Real, PCR RÁPIDO y extra resistente
SP - ACCESORIOS OPCIONALES: la tira de tapas
SP - APROPIADO PARA: PCR
SP - APROPIADO PARA: PCR Real Time
SP - APROPIADO PARA: PCR Real Time FAST
SP - APROPIADO PARA: almacenamiento
SP - PAREDES: delgado 0,30 mm
SP - TRATAMIENTO DE LA SEÑAL: las paredes de cristal esmerilado
SP - INTENSIDAD DE SEÑAL: 125
SP - A LA POSIBLE INTERFERENCIA DE LA SEÑAL: 70
SP - RESISTENTE: la congelación a -20°C
SP - EVAPORACIÓN DE LA MUESTRA: el menor de 3%
SP - UNIÓN: Menor del ADN 0.3%
SP - UNIÓN: las proteínas de menos de 0,8%
SP - COMPATIBLE CON: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE CON: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml



SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-Opticon System-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml