



Strip 8x0.1ml (LP) extra-rigide - bianche

Cf 120pz - Strip 8x0.1ml (LP) extra-rigide - bianche. Cf 120pz

Product codes:

Reference: 12 954 9671



Product gallery:



Product description:

Strip da 8 microtubi compatibile con molti strumenti che accettano plastica per biologia molecolare di tipo FAST - Low Profile

- Disegnata per applicazioni in PCR e PCR Real Time , compatibile anche con apparecchi come ABI / Life Technologies, Roche, Bio-Rad, Eppendorf e altri
- Striscia robusta non flessibile
- Versione Low Profile da 0.1ml
- Colore bianco
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

Product features:

- SP - TIPO DE: strip da 8 microtubi
- SP - CAP: senza
- SP - VOLUMEN: 0.1ml
- SP - CAPACIDAD: 0.2ml
- SP - POZOS: 8
- SP - PERFIL: tipo de baja RÁPIDO/0.1



SP - ALTURA: 15.6 / 16.6mm ca
SP - FORMA: tira resistente
SP - FONDO: cónico
SP - MATERIAL: PP
SP - APLICACIONES: PCR
SP - APLICACIONES: PCR Real Time
SP - APLICACIONES: PCR Real Time FAST
SP - ESTERILIZABLE EN AUTOCLAVE: uested
SP - COLOR: bianco
SP - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - ESTERILIDAD: no
SP - CONFORMIDAD: no sujeto a CE
SP - TEMPERATURA MAX: 100°C
SP - TEMPERATURA MÍNIMA: -30°C
EL SP - PAQUETE: sobres
SP - CARACTERÍSTICAS: strip da 8x0.1 per PCR Real Time FAST e PCR a striscia robusta - bianche
SP - ACCESORIOS OPCIONALES: griglia assemblaggio bianca
SP - ACCESORIOS OPCIONALES: griglia assemblaggio gialla
SP - APROPIADO PARA: PCR
SP - APROPIADO PARA: PCR Real Time
SP - APROPIADO PARA: PCR Real Time FAST
SP - PAREDES: delgado 0,30 mm
SP - TRATAMIENTO DE LA SEÑAL: pareti bianche
SP - INTENSIDAD DE SEÑAL: 180
SP - A LA POSIBLE INTERFERENCIA DE LA SEÑAL: inferiore a 2
SP - RESISTENTE: la congelación a -20°C
SP - EVAPORACIÓN DE LA MUESTRA: el menor de 3%
SP - UNIÓN: Menor del ADN 0.3%
SP - UNIÓN: las proteínas de menos de 0,8%
SP - COMPATIBLE CON: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE CON: Applied Biosystems-9800 Fast Thermal Cyclers-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml



SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-Opticon System-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE CON: Roche-LightCycler® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml



SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Roche-LightCycler® 480 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
SP - COMPATIBLE CON: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml