



Strip 8x0.2ml (RP) extra-rigide - bianche

Cf 120pz - Strip 8x0.2ml (RP) extra-rigide - bianche. Cf 120pz



Product codes:

Reference: 12 954 7791

Product gallery:

Product description:

Strip da 8 microtubi compatibile con molti strumenti per PCR e PCR Real Time che accettano plastica per biologia molecolare a profilo regolare e alcuni modelli di sequenziatori

- Disegnata per applicazioni in PCR e PCR Real Time , compatibile anche con apparecchi come ABI / Life Technologies, Bio-Rad, Eppendorf e altri
- Striscia robusta non flessibile
- Versione profilo regolare da 0.2ml
- Colore bianco
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

Product features:

SP - TIPO DE: strip da 8 microtubi



- SP - CAP: senza
- SP - VOLUMEN: 0.2ml
- SP - CAPACIDAD: 0.2ml
- SP - POZOS: 8
- SP - PERFIL: ajustar
- SP - ALTURA: 20.9 / 21.8mm ca
- SP - FORMA: tira resistente
- SP - FONDO: cónico
- SP - MATERIAL: PP
- SP - APLICACIONES: PCR
- SP - APLICACIONES: PCR Real Time
- SP - ESTERILIZABLE EN AUTOCLAVE: uested
- SP - COLOR: bianco
- SP - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
- SP - ESTERILIDAD: no
- SP - CONFORMIDAD: no sujeto a CE
- SP - TEMPERATURA MAX: 100°C
- SP - TEMPERATURA MÍNIMA: -30°C
- EL SP - PAQUETE: sobres
- SP - CARACTERÍSTICAS: strip da 8x0.1 per PCR Real Time e PCR a striscia robusta - bianche
- SP - ACCESORIOS OPCIONALES: la cuadrícula de la asamblea azul
- SP - APROPIADO PARA: PCR
- SP - APROPIADO PARA: PCR Real Time
- SP - PAREDES: delgado 0,30 mm
- SP - TRATAMIENTO DE LA SEÑAL: pareti bianche
- SP - INTENSIDAD DE SEÑAL: 180
- SP - A LA POSIBLE INTERFERENCIA DE LA SEÑAL: inferiore a 2
- SP - RESISTENTE: la congelación a -20°C
- SP - EVAPORACIÓN DE LA MUESTRA: el menor de 3%
- SP - UNIÓN: Menor del ADN 0.3%
- SP - UNIÓN: las proteínas de menos de 0,8%
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
- SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
- SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
- SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
- SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
- SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
- SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
- SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
- SP - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml



- SP - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
- SP - COMPATIBLE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
- SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
- SP - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-310-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-3100-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-3130-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-3500 Genetic Analyzer-8-capillary
- SP - COMPATIBLE CON: Applied Biosystems-3700-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-3730-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
- SP - COMPATIBLE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
- SP - COMPATIBLE CON: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
- SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml



SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-SimpliAmp™-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-Solid3-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
SP - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
SP - COMPATIBLE CON: Stratagene-Mx3000P™-96 x 0.2ml
SP - COMPATIBLE CON: Stratagene-Mx3005P™-96 x 0.2ml
SP - COMPATIBLE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TOptical-96 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml



-
- SP - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
 - SP - COMPATIBLE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
 - SP - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
 - SP - COMPATIBLE CON: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
 - SP - COMPATIBLE CON: Applied Biosystems-3500xL Genetic Analyzer-24-capillary
 - SP - COMPATIBLE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
 - SP - COMPATIBLE CON: Applied Biosystems-MiniAmp Plus Thermal Cycler-96 x 0.2ml
 - SP - COMPATIBLE CON: Applied Biosystems-MiniAmp Thermal Cycler-96 x 0.2ml
 - SP - COMPATIBLE CON: Stratagene-Mx4000™-96 x 0.2ml