



1.88 inch wide Crystal Clear Sealing Tape, with cutter - 1.88 inch x 43.7 yard roll, with cutter - 1/ea

Sealing tape used to seal sitting drop crystallization experiments

Product codes:

Reference: HR 3511

Product gallery:



Plate	96 Wells	48 Wells	96 Wells
Corning Costar	Yes	Yes	Yes
Hampton Research CrystalClear	Yes	Yes	Yes
Hampton Research CrystalClear 96	Yes	Yes	Yes
Hampton Research CrystalClear 48	Yes	Yes	Yes
Hampton Research CrystalClear 96	Yes	Yes	Yes
Hampton Research CrystalClear 48	Yes	Yes	Yes
Hampton Research CrystalClear 96	Yes	Yes	Yes
Hampton Research CrystalClear 48	Yes	Yes	Yes
Hampton Research CrystalClear 96	Yes	Yes	Yes
Hampton Research CrystalClear 48	Yes	Yes	Yes

Product description:

- Optically clear
- Compatible with a wide range of crystallization reagents
- UV compatible
- Crystal Clear Mini Sealing Tape reseals and patches 96 well sitting drop plates
- Crystal Clear Mini Sealing Tape seals Cryschem S and Cryschem M plate one row at a time

These optically clear tapes are compatible with protein crystallization reagents. The polypropylene Crystal Clear Sealing Tapes use acrylic solvent based adhesives which are compatible with aqueous crystallization reagents.

Catalog number HR3-511 is a 1.88 inch (48 mm) wide, 3 mil tape on a 43.7 yard (40 M) roll with a 1.5 inch core and is supplied with a green dispenser / cutter.

Catalog number HR4-511 is a 1.88 inch wide, 3 mil tape on a 109.3 yard roll with a 3 inch core and no dispenser / cutter.

Two strips of the HR3-511 or HR4-511 will seal a Corning, Cryschem, CrystalClear Strip, Greiner, Intelliplate, Linbro, Swissci (MRC), VDX and VDXm plate.

Catalog number HR4-506 is a 3 inch wide, 3 mil tape on a 54.86 yard roll with a 3 inch core and no dispenser / cutter.



One strip of the HR4-506 will seal a Corning, Cryschem M, CrystalClear Strip, Greiner, Intelliplate, Swissci (MRC) and VDXm plate.

Crystal Clear Mini Sealing Tape HR4-508 is a 0.75 inch (19 mm) wide, 3 mil tape on a 650 inch (16.5 M) roll and is supplied with a clear dispenser / cutter. The tape is useful for resealing or patching wells after mounting crystals. The tape is also wide enough to seal a single row of a Cryschem S or Cryschem M plate.

See chart for tape and plate compatibility.

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

Product features:

CRF - TIPO: Crystal Clear Sealing Tape;