



CrystalWand Magnetic, Straight - 1/ea

CrystalCap Magnetic & CrystalCap Copper Magnetic handling tool



Codes produits :

Référence HR 4729

Galerie de produits :



Description du produit :

The CrystalWand Magnetic is designed to be used exclusively with the CrystalCap, CrystalCap ALS, and CrystalCap SPINE systems during transfer of the caps from the vial to the goniometer head and from the goniometer to the vial. The 6 1/2" (165 mm) long chrome plated steel wand features a plastic housing enclosing a spring tensioned plunger that when depressed, moves a non-magnetic steel platform away from the magnet housed in the end of the wand. This causes the steel CrystalCap to detach from the wand and be placed readily into the vial. When the platform is retracted, the wand securely holds the steel CrystalCap Magnetic.

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

caractéristiques du produit :

CRF - TIPO: CrystalWand Magnetic