



Strips 8 tappi ottici piani anti-impronta

Cf. 120 strips da 8 tappi - Strips di 8 tappi ottici piani anti-impronta

Codes produits :

Référence 12 864 8451



Galerie de produits :



Description du produit :

Strips di 8 tappi ottici piani anti-impronta

La superficie rientrante evita che l'area ottica venga "contaminata" da impronte accidentali

Utilizzabili in tutte le procedure in PCR Real Time

Pareti sottili ed uniformi, ottimi per chiudere le piastre e le strips per le quali sono indicate come accessorio.

Ottimizzate per PCR Real Time

Cf. 120 strips da 8 tappi



caractéristiques du produit :

FLM - TYPE: la bande de casquettes
FLM - COMPATIBILITÉ: strip PCR
FLM - COMPATIBILITÉ: plaques PCR
FLM - COMPATIBILITÉ: plaques de PCR en Temps Réel
FLM - COMPATIBILITÉ: strip PCR Real-Time
FLM - MATÉRIEL: PP
FLM - PUIITS: 8
FLM - PURETÉ: DNA / RNA / DNase / Rnase / Pyrogen FREE
FLM - STÉRILITÉ: pas de
FLM - FORME: réentrant
FLM - FORME: bande robuste casquettes
FLM - TEMPÉRATURE MAX: 100°C
FLM - TEMPÉRATURE MINIMALE: -30°C
FLM - OPTIQUE DE LA ZONE DE: 12.6mm
FLM - DOMAINE ACCESSIBLE EN ÉCRITURE: vous
FLM - AUTOCLAVABLE: vous
FLM - CODAGE COULEUR: pas de
FLM - ÉTIQUETTE: pas de
FLM - CONFORMITÉ: pas soumis à CE
FLM - CARACTÉRISTIQUES: strip da 8 tappi piani anti-interferenza extra-robusti
FLM - ADAPTÉ POUR: la biologie moléculaire
FLM - ADAPTÉ POUR: PCR
FLM - ADAPTÉ POUR: PCR Real Time
FLM - ADAPTÉ POUR: PCR Real Time FAST
FLM - ADAPTÉ POUR: de stockage
FLM - ADAPTÉ POUR: le stockage à -20°C
FLM - COMPATIBLE AVEC: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBLE AVEC: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBLE AVEC: Applied Biosystems-9800 Fast Thermal Cyclers-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml



FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
FLM - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-SimpliAmp™-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-Solid3-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
FLM - COMPATIBLE AVEC: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBLE AVEC: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBLE AVEC: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Deep Well
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Well Fast
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-iCycler-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Deep Well
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Well Fast
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-T100-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler-96-77 96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler gradient-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml



FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler pro-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Mastercycler pro S-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Thermomixer C-96 x 0.2ml
FLM - COMPATIBLE AVEC: Eppendorf-Thermomixer comfort-96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-Opticon System-96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBLE AVEC: Roche-cobas z 480-96 x 0.2ml
FLM - COMPATIBLE AVEC: Roche-LightCycler ® 480 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBLE AVEC: Stratagene-Mx3000P™-96 x 0.2ml
FLM - COMPATIBLE AVEC: Stratagene-Mx3005P™-96 x 0.2ml
FLM - COMPATIBLE AVEC: Stratagene-Mx4000™-96 x 0.2ml
FLM - COMPATIBLE AVEC: Stratagene-Robocycler Gradient 96-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TGradient-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TOptical-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TProfessional-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TRobot-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-UNO II-96 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-iCycler-48 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-Dual 48/48 Fast



FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: Roche-COBAS Taqman 48 Analyzer-24 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TPersonal-48 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
FLM - COMPATIBLE AVEC: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
FLM - COMPATIBLE AVEC: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
FLM - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBLE AVEC: Roche-LightCycler® Nano-32 x 0.2ml
FLM - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler-48 x 0.2ml