



Microtubi qPCR da 0,2ml


Cf. 1000 pezzi - Microtubi da 0,2ml per qPCR, profilo regolare



Codes produits :

Référence 12 879 7601

Galerie de produits :



Description du produit :

Microtubi da 0,2ml, profilo regolare, eccellenti per PCR Real-Time

Compatibili con molti apparecchi per PCR Real-Time che accettano provettine da 0,2ml a profilo regolare.

Non compatibile con apparecchi FAST.

Cf. 1000pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Leggenda:

Grado qPCR 1: sviluppato per segnale ottimo in PCR Real-Time

Grado qPCR 2: sviluppato per avere il segnale al top, esente da interferenze, in PCR Real-Time



caractéristiques du produit :

VTT - CAP: plat
MTB - CONFORMITÉ: pas soumis à CE
VTT - TYPE: PCR
MTB - VOLUME: 0.2ml
MTB - CAPACITÉ: 0.2ml
VTT - CAP FORABLE: pas de
MTB - INSÉRÉ CAP: pas de
VTT DE LA RÉGION OPTIQUE: 12.6mm
MTB - FORME: ajuster
MTB - FONDS: conique
MTB - AUTONOME: pas de
VTT - MATÉRIEL: PP
VTT TAILLE: pas de
MTB - OPACITÉ: légèrement opaque
VTT - ZONE INSCRIPTIBLE: pas de
MTB - PLAGE DE TEMPÉRATURE: -30°C à +100°C
MTB - TEMPÉRATURE MAX: 100°C
VTT DE LA PURETÉ: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - AUTOCLAVABLE: vous
VTT - LA STÉRILITÉ: pas de
VTT - TAILLE: h 20.9/21.8mm
VTT - CARACTÉRISTIQUES: microtubes de 0,2 ml avec bouchon plat pour la PCR et PCR en Temps Réel
MTB - ADAPTÉ POUR: la préparation de l'échantillon
MTB - ADAPTÉ POUR: de stockage
MTB - ADAPTÉ POUR: PCR
MTB - ADAPTÉ POUR: PCR Real Time
MTB - RÉSISTANT: la congélation à -20°C
MTB - EMBALLAGE: enveloppe en vrac
MTB - MURS: mince de 0,3 mm
MTB - L'ÉVAPORATION DE L'ÉCHANTILLON: le moindre de 3%
VTT DE LIAISON: ADN mineur 0.3%
VTT DE LIAISON: les protéines de moins de 0,8%
MTB - INDICATEUR D'INTENSITÉ DU SIGNAL: 125
VTT - LES INTERFÉRENCES DE SIGNAL: 70
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Deep Well
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Well Fast
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml



MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-iCycler-48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-iCycler-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Deep Well
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Well Fast
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-T100-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler-96-77 96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler gradient-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler pro-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler pro S-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Thermomixer C-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Thermomixer comfort-96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml



MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TGradient-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TPersonal-48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TRobot-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-UNO II-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-310-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-3100-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-3130-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-3700-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-3730-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBLE AVEC: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBLE AVEC: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-SimpliAmp™-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml



MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBLE AVEC: Roche-COBAS Taqman 48 Analyzer-24 x 0.2ml
MTB - COMPATIBLE AVEC: Stratagene-Mx3000P™-96 x 0.2ml
MTB - COMPATIBLE AVEC: Stratagene-Mx3005P™-96 x 0.2ml
MTB - COMPATIBLE AVEC: Stratagene-Mx4000™-96 x 0.2ml
MTB - COMPATIBLE AVEC: Stratagene-Robocycler Gradient 96-96 x 0.2ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TOptical-96 x 0.2ml