




Microtubi PCR da 0,5ml tappo piatto, graduati

Cf. 1000 pezzi - Microtubi da 0,5ml per PCR, tappo piatto, graduati



Codes produits :

Référence 12 879 3051

Galerie de produits :



Description du produit :

Microtubi da 0,5ml, tappo piatto, graduati.

Compatibili con molti apparecchi per PCR che accettano provettine da 0,5ml con tappo piatto.

Profilo regolare.

Cf. 1000pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

caractéristiques du produit :

VTT - CAP: plat



MTB - CONFORMITÉ: pas soumis à CE
VTT - TYPE: PCR
VTT - TYPE: Eppendorf
MTB - VOLUME: 0.5ml
MTB - CAPACITÉ: 0.5ml
VTT - CAP FORABLE: pas de
MTB - INSÉRÉ CAP: pas de
VTT DE LA RÉGION OPTIQUE: 5mm
MTB - FORME: ajuster
MTB - FONDS: conique
MTB - AUTONOME: pas de
VTT - MATÉRIEL: PP
VTT TAILLE: vous
MTB - OPACITÉ: ultra-transparent
VTT - ZONE INSCRIPTIBLE: pas de
MTB - PLAGE DE TEMPÉRATURE: -30°C a +100°C
MTB - TEMPÉRATURE MAX: 100°C
VTT DE LA PURETÉ: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - AUTOCLAVABLE: vous
VTT - LA STÉRILITÉ: pas de
VTT - TAILLE: h 30mm ca
VTT - CARACTÉRISTIQUES: microtubi da 0.5ml con tappo piatto graduati per PCR e PCR Real-Time
MTB - ADAPTÉ POUR: la préparation de l'échantillon
MTB - ADAPTÉ POUR: de stockage
MTB - ADAPTÉ POUR: PCR
MTB - ADAPTÉ POUR: PCR Real Time
MTB - RÉSISTANT: la congélation à -20°C
MTB - EMBALLAGE: enveloppe en vrac
MTB - MURS: mince de 0,3 mm
VTT DE LIAISON: les protéines de moins de 0,8%
VTT DE LIAISON: ADN mineur 0.3%
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler-48 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler Combi-96-77 77 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler-20 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler-20 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TGradient-48 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TPersonal-20 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TPersonal Combi-48-18 18 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional-60 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional Standard-60 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO-30 x 0.5ml



MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-18 18 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TRIO Thermoblock-20 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-UNO II-48 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-UNO-thermoblock-40 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-GeneCycler-12 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 480-48 x 0.5ml
MTB - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9700-60 x 0.5ml
MTB - COMPATIBLE AVEC: Applied Biosystems-Veriti™ 60-Well Thermal Cycler-60 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-60 Alpha unit 60 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-iCycler-60 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-MJ Mini-48-12 12 x 0.5ml
MTB - COMPATIBLE AVEC: Bio-Rad Laboratories-PTC-100-60 x 0.5ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler-96-77 77 x 0.5ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler gradient-96-77 77 x 0.5ml
MTB - COMPATIBLE AVEC: Eppendorf-Mastercycler Personal-25-16 16 x 0.5ml
MTB - COMPATIBLE AVEC: Eppendorf-Thermomixer C-24 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-60 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-150 MiniCycler-16 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-30-30 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 30 x 0.5ml
MTB - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-60 Alpha unit 60 x 0.5ml