



Strips 8 pozzetti qPCR da 0,1ml

Cf. 120 pezzi - Microtubi da 0,1ml per qPCR, profilo basso

Codes produits :

Référence 12 954 8391



Galerie de produits :



Description du produit :

Strips da 0,1ml, profilo basso, eccellenti per PCR Real-Time.

La scelta migliore per gli apparecchi ABI e BioRad per PCR Real-Time a basso profilo.

Tutta la serie di strips possiede una forma asimmetrica come riferimento orientativo, un'area scrivibile e dei fori per fissarli sulla griglia atta ad ospitare le strips per creare la propria piastra.

Cf. 120pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

Leggenda:

Grado qPCR 1: sviluppato per segnale ottimo in PCR Real-Time

Grado qPCR 2: sviluppato per avere il segnale al top, esente da interferenze, in PCR Real-Time



caractéristiques du produit :

SP - TYPE: strip da 8 microtubi

SP - CAP: senza

SP - VOLUME: 0.1ml

SP - CAPACITÉ: 0.2ml

SP - PUIITS: 8

SP - PROFIL: bas de type RAPIDE/0.1

SP - HAUTEUR: 15.6 / 16.6mm ca

SP - FORME: bande robuste

SP - FONDS: conique

SP - MATÉRIEL: PP

SP - APPLICATIONS: PCR

SP - APPLICATIONS: PCR Real Time

SP - APPLICATIONS: PCR Real Time FAST

SP - AUTOCLAVABLE: vous

SP - COULEUR: bianco

SP - PURETÉ: DNA / RNA / DNase / RNase / Pyrogen FREE

SP - STÉRILITÉ: pas de

SP - CONFORMITÉ: pas soumis à CE

SP - TEMPÉRATURE MAX: 100°

SP - TEMPÉRATURE MINIMALE: -30°C

LE SP - PAQUET: enveloppe

SP - CARACTÉRISTIQUES: strip da 8x0.1 per PCR Real Time FAST e PCR extra robuste bianco

SP - ACCESSOIRES EN OPTION: la bande de casquettes

SP - ACCESSOIRES EN OPTION: griglia assemblaggio bianca

SP - ADAPTÉ POUR: PCR

SP - ADAPTÉ POUR: PCR Real Time

SP - ADAPTÉ POUR: PCR Real Time FAST

SP - ADAPTÉ POUR: de stockage

SP - MURS: mince 0,30 mm

SP - TRAITEMENT DU SIGNAL: parietal blanche

SP - INDICATEUR D'INTENSITÉ DU SIGNAL: 180

SP - POSSIBLE LES INTERFÉRENCES DE SIGNAL: inferiore a 2

SP - RÉSISTANT: la congélation à -20°C

SP - ÉVAPORATION DE L'ÉCHANTILLON: le moindre de 3%

SP - LIAISON: ADN mineur 0.3%

SP - LIAISON: les protéines de moins de 0,8%

SP - COMPATIBLE AVEC: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml

SP - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml

SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast

SP - COMPATIBLE AVEC: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml

SP - COMPATIBLE AVEC: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml

SP - COMPATIBLE AVEC: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml

SP - COMPATIBLE AVEC: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml



SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-Opticon System-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Roche-LightCycler® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic-96 x 0.2ml



SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Roche-LightCycler® 480 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml