



Strip 8x0.1ml (LP) extra-rigide - blanche

Cf 120pz - Strip 8x0.1ml (LP) extra-rigide - blanche. Cf 120pz

Codes produits :

Référence 12 954 9671



Galerie de produits :



Description du produit :

Strip da 8 microtubi compatibile con molti strumenti che accettano plastica per biologia molecolare di tipo FAST - Low Profile

- Disegnata per applicazioni in PCR e PCR Real Time , compatibile anche con apparecchi come ABI / Life Technologies, Roche, Bio-Rad, Eppendorf e altri
- Striscia robusta non flessibile
- Versione Low Profile da 0.1ml
- Colore bianco
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

caractéristiques du produit :

SP - TYPE: strip da 8 microtubi

SP - CAP: senza

SP - VOLUME: 0.1ml

SP - CAPACITÉ: 0.2ml

SP - PUIITS: 8

SP - PROFIL: bas de type RAPIDE/0.1



SP - HAUTEUR: 15.6 / 16.6mm ca
SP - FORME: bande robuste
SP - FONDS: conique
SP - MATÉRIEL: PP
SP - APPLICATIONS: PCR
SP - APPLICATIONS: PCR Real Time
SP - APPLICATIONS: PCR Real Time FAST
SP - AUTOCLAVABLE: vous
SP - COULEUR: bianco
SP - PURETÉ: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STÉRILITÉ: pas de
SP - CONFORMITÉ: pas soumis à CE
SP - TEMPÉRATURE MAX: 100°C
SP - TEMPÉRATURE MINIMALE: -30°C
LE SP - PAQUET: enveloppe
SP - CARACTÉRISTIQUES: strip da 8x0.1 per PCR Real Time FAST e PCR a striscia robusta - bianche
SP - ACCESSOIRES EN OPTION: griglia assemblaggio bianca
SP - ACCESSOIRES EN OPTION: griglia assemblaggio gialla
SP - ADAPTÉ POUR: PCR
SP - ADAPTÉ POUR: PCR Real Time
SP - ADAPTÉ POUR: PCR Real Time FAST
SP - MURS: mince 0,30 mm
SP - TRAITEMENT DU SIGNAL: parietali bianche
SP - INDICATEUR D'INTENSITÉ DU SIGNAL: 180
SP - POSSIBLE LES INTERFÉRENCES DE SIGNAL: inferiore a 2
SP - RÉSISTANT: la congélation à -20°C
SP - ÉVAPORATION DE L'ÉCHANTILLON: le moindre de 3%
SP - LIAISON: ADN mineur 0.3%
SP - LIAISON: les protéines de moins de 0,8%
SP - COMPATIBLE AVEC: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE AVEC: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBLE AVEC: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
SP - COMPATIBLE AVEC: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler-96-77 96 x 0.2ml



SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBLE AVEC: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-Opticon System-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBLE AVEC: Roche-LightCycler® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBLE AVEC: Stratagene-Robocycler Gradient 96-96 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBLE AVEC: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml



SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBLE AVEC: Roche-LightCycler® 480 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBLE AVEC: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
SP - COMPATIBLE AVEC: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
SP - COMPATIBLE AVEC: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
SP - COMPATIBLE AVEC: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
SP - COMPATIBLE AVEC: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE AVEC: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBLE AVEC: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml