

Strip 8 tappi ottici piani anti-impronta extra robusti

Cf. 120 strips da 8 tappi - Strips di 8 tappi ottici piani anti-impronta

Product codes:

Reference: 12 864 8451



Product gallery:



Product description:

Strips di 8 tappi ottici piani anti-impronta

La superficie rientrante evita che l'area ottica venga "contaminata" da impronte accidentali

Utilizzabili in tutte le procedure in PCR Real Time

Pareti sottili ed uniformi, ottimi per chiudere le piastre e le strips per le quali sono indicate come accessorio.

Ottimizzate per PCR Real Time

Cf. 120 strips da 8 tappi

COMPATIBILI CON:

Applied Biosystems 7300 Real-Time PCR System 96 x 0.2ml Applied Biosystems 7500 Fast Real-Time PCR System 96 x 0.2ml Fast Applied Biosystems 7500 Real-Time PCR System 96 x 0.2ml Standard Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Fast Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Standard Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml Applied Biosystems ABI PRISM 7000 SDS 96 x 0.2ml Applied Biosystems ABI PRISM 7700 SDS 96 x 0.2ml Applied Biosystems GeneAmp PCR System 2400 24 x 0.2ml Applied Biosystems GeneAmp PCR System 2700 96 x 0.2ml Applied Biosystems GeneAmp PCR System 2720 96 x 0.2ml Applied Biosystems GeneAmp PCR System 5700 96 x 0.2ml Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml Applied Biosystems GeneAmp PCR System 9700 96 x 0.2ml Applied Biosystems GeneAmp PCR System 9700 Dual 96-Well 96 x 0.2ml Applied Biosystems GeneAmp PCR System 9700 Aluminum 96 x 0.2ml Applied Biosystems GeneAmp PCR System 9700 Gold 96 x 0.2ml Applied Biosystems ProFlex™ PCR System 32 x 0.2ml Applied Biosystems ProFlex™ PCR System 96 x 0.2ml Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 12k Flex 96 x 0.2ml Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.2ml Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.2ml Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 6 Flex 96 x 0.2ml Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 7 Flex 96 x 0.2ml Applied Biosystems SimpliAmp™ 96 x 0.2ml Applied Biosystems Solid3 96 x 0.2ml Applied Biosystems StepOnePlus™ Real-Time PCR System 96 x 0.1ml Applied Biosystems StepOne™ Real-Time PCR System 48 x 0.1ml Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml Applied Biosystems Veriti™ 96-Well Thermal Cycler 96 x 0.2ml Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml Applied Biosystems Veriti™ Dx 96-Well Thermal Cycler 96 x 0.2ml Applied Biosystems ViiA™ 7 Real-Time PCR System 96 x 0.2ml Applied Biosystems ViiA™ 7 Real-Time PCR System Fast 96 x 0.1ml Bio-Rad Laboratories C1000 96-Well Fast Bio-Rad Laboratories C1000 Dual 48/48 Fast Bio-Rad Laboratories C1000 Touch 96-Deep Well Bio-Rad Laboratories C1000 Touch 96-Well Fast Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml Bio-Rad Laboratories DNA Engine 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Disciple 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml Bio-Rad Laboratories DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml Bio-Rad Laboratories GeneCycler 24 x 0.2ml Bio-Rad Laboratories iCycler 48 x 0.2ml Bio-Rad Laboratories iCycler 96 x 0.2ml Bio-Rad Laboratories iQ Real-Time PCR 96 x 0.2ml Bio-Rad Laboratories iQ5 Real-Time PCR 96 x 0.2ml Bio-Rad

Laboratories MiniOpticon 48-12 48 x 0.2ml Bio-Rad Laboratories MJ Mini 48-12 48 x 0.2ml Bio-Rad Laboratories MyCycler 96 x 0.2ml Bio-Rad Laboratories MyiQ Real-Time PCR 96 x 0.2ml Bio-Rad Laboratories MyiQ2 Real-Time PCR 96 x 0.2ml Bio-Rad Laboratories Opticon 2 96 x 0.2ml Bio-Rad Laboratories PTC-100 96 x 0.2ml Bio-Rad Laboratories S1000 96-Deep Well Bio-Rad Laboratories S1000 96-Well Fast Bio-Rad Laboratories S1000 Dual 48/48 Fast Bio-Rad Laboratories T100 96 x 0.2ml Eppendorf Mastercycler 96 x 0.2ml Eppendorf Mastercycler 96-77 96 x 0.2ml Eppendorf Mastercycler ep gradient 96 x 0.2ml Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml Eppendorf Mastercycler ep realplex 96 x 0.2ml Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml Eppendorf Mastercycler gradient 96 x 0.2ml Eppendorf Mastercycler gradient 96-77 96 x 0.2ml Eppendorf Mastercycler nexus 96 x 0.2ml Eppendorf Mastercycler nexus eco 96 x 0.2ml Eppendorf Mastercycler nexus gradient 96 x 0.2ml Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml Eppendorf Mastercycler nexus GSX1 96 x 0.2ml Eppendorf Mastercycler nexus GSX1e 96 x 0.2ml Eppendorf Mastercycler nexus SX1 96 x 0.2ml Eppendorf Mastercycler nexus SX1e 96 x 0.2ml Eppendorf Mastercycler pro 96 x 0.2ml Eppendorf Mastercycler pro S 96 x 0.2ml Eppendorf Thermomixer C 96 x 0.2ml Eppendorf Thermomixer comfort 96 x 0.2ml MJ Research Inc. Opticon 2 System 96 x 0.2ml MJ Research Inc. Opticon System 96 x 0.2ml MJ Research Inc. PTC-100 96 x 0.2ml MJ Research Inc. PTC-100 96 x 0.2ml AgV MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-200 DNA Engine 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml Roche COBAS Taqman 48 Analyzer 24 x 0.2ml Roche cobas z 480 96 x 0.2ml Roche LightCycler® 480 Real-Time PCR System 96 x 0.2ml Roche LightCycler® 96 Real-Time PCR System 96 x 0.2ml Roche LightCycler® Nano 32 x 0.2ml Stratagene Mx3000P™ 96 x 0.2ml Stratagene Mx3005P™ 96 x 0.2ml Stratagene Mx4000™ 96 x 0.2ml Stratagene Robocycler Gradient 96 96 x 0.2ml Whatman Biometra T1 Thermocycler 96 x 0.2ml Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml Whatman Biometra T3 Thermocycler 48 x 0.2ml Whatman Biometra T3 Thermocycler Combi 48-18 48 x 0.2ml Whatman Biometra T3000 Thermocycler 48 x 0.2ml Whatman Biometra T3000 Thermocycler Combi 48-18 48 x 0.2ml Whatman Biometra TGradient 96 x 0.2ml Whatman Biometra TOptical 96 x 0.2ml Whatman Biometra TPersonal 48 x 0.2ml Whatman Biometra TPersonal Combi 48-18 48 x 0.2ml

Whatman Biometra TProfessional 96 x 0.2ml Whatman Biometra TProfessional Basic 96 x 0.2ml
Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml Whatman Biometra TProfessional
Standard Gradient 96 x 0.2ml Whatman Biometra TProfessional TRIO 48 x 0.2ml Whatman
Biometra TProfessional TRIO combi 48-30 48 x 0.2ml Whatman Biometra TRobot 96 x 0.2ml
Whatman Biometra UNO II 96 x 0.2ml Whatman Biometra UNO-thermoblock 96 x 0.2ml

Product features:

FLM - TIPO: strip di tappi
FLM - COMPATIBILITA': strip PCR
FLM - COMPATIBILITA': piastre PCR
FLM - COMPATIBILITA': piastre PCR Real-Time
FLM - COMPATIBILITA': strip PCR Real-Time
FLM - MATERIALE: PP
FLM - POZZETTI: 8
FLM - PUREZZA: DNA / RNA / DNase / Rnase / Pyrogen FREE
FLM - STERILITA': no
FLM - FORMA: rientrante
FLM - FORMA: striscia robusta di tappi
FLM - TEMPERATURA MAX: 100°C
FLM - TEMPERATURA MINIMA: -30°C
FLM - AREA OTTICA: 12.6mm
FLM - AREA SCRIVIBILE: si
FLM - AUTOCLAVABILE: si
FLM - CODIFICA COLORE: no
FLM - ETICHETTA: no
FLM - CONFORMITA': non soggetto a CE
FLM - PECULIARITA': strip da 8 tappi piani anti-interferenza extra-robusti
FLM - IDONEE PER: biologia molecolare
FLM - IDONEE PER: PCR
FLM - IDONEE PER: PCR Real Time
FLM - IDONEE PER: PCR Real Time FAST
FLM - IDONEE PER: stoccaggio
FLM - IDONEE PER: stoccaggio a -20°C
FLM - COMPATIBILE CON: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBILE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
FLM - COMPATIBILE CON: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml

FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.1ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 12k Flex-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.1ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.1ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 6 Flex-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.1ml Fast
FLM - COMPATIBILE CON: Applied Biosystems-QuantStudio™ 7 Flex-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-SimpliAmp™-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-Solid3-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-StepOnePlus™ Real-Time PCR System-96 x 0.1ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ViiA™ 7 Real-Time PCR System Fast-96 x 0.1ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Well Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Well Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml

FLM - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Thermomixer C-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-Opticon System-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
FLM - COMPATIBILE CON: Roche-cobas z 480-96 x 0.2ml
FLM - COMPATIBILE CON: Roche-LightCycler ® 480 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Mx3000P™-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Mx3005P™-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Mx4000™-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TGradient-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TOptical-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TRobot-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-UNO II-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml

FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: Roche-COBAS Taqman 48 Analyzer-24 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TPersonal-48 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-StepOne™ Real-Time PCR System-48 x 0.1ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
FLM - COMPATIBILE CON: Roche-LightCycler® Nano-32 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml