

Strip 12 tappi ottici piani anti-impronta extra robusti

Strips di 12 tappi ottici piani anti-impronta


Cf. 80 strip



Product codes:

Reference: 12 864 4011

Product gallery:



Indented qPCR Cap
the "optical area" is
indented to the cap
surface which
prevents "touching"
the optical area

Product description:

Strips di 12 tappi ottici piani anti-impronta

La superficie rientrante evita che l'area ottica venga "contaminata" da impronte accidentali

Utilizzabili in tutte le procedure in PCR

Pareti sottili ed uniformi, ottimi per chiudere le piastre e le strips per le quali sono indicate come accessorio.

Ottimizzate per PCR

Cf. 80 strip da 12 tappi

IDEONEE PER UTILIZZI CON I SEGUENTI APPARECCHI.

MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml
Stratagene	Mx3000P™	96 x 0.2ml
Stratagene	Mx3005P™	96 x 0.2ml
Stratagene	Mx4000™	96 x 0.2ml
Stratagene	Robocycler Gradient 96	96 x 0.2ml
Whatman Biometra	T1 Thermocycler	96 x 0.2ml
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml
Whatman Biometra	TGradient	96 x 0.2ml
Whatman Biometra	TOptical	96 x 0.2ml
Whatman Biometra	TProfessional	96 x 0.2ml
Whatman Biometra	TProfessional Basic	96 x 0.2ml
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml
Whatman Biometra	TProfessional Standard Gradient	96 x 0.2ml
Whatman Biometra	TRobot	96 x 0.2ml
Whatman Biometra	UNO II	96 x 0.2ml
Whatman Biometra	UNO-thermoblock	96 x 0.2ml

Product features:

- FLM - TIPO: strip di tappi
- FLM - COMPATIBILITA': strip PCR
- FLM - COMPATIBILITA': piastre PCR
- FLM - COMPATIBILITA': piastre PCR Real-Time
- FLM - COMPATIBILITA': strip PCR Real-Time
- FLM - MATERIALE: PP
- FLM - POZZETTI: 12
- FLM - PUREZZA: DNA / RNA / DNase / Rnase / Pyrogen FREE
- FLM - STERILITA': no
- FLM - FORMA: rientrante
- FLM - FORMA: striscia robusta di tappi
- FLM - TEMPERATURA MAX: 100°C
- FLM - TEMPERATURA MINIMA: -30°C
- FLM - AREA OTTICA: 12.6mm
- FLM - AREA SCRIVIBILE: si
- FLM - AUTOCLAVABILE: si
- FLM - CODIFICA COLORE: no
- FLM - ETICHETTA: no
- FLM - CONFORMITA': non soggetto a CE
- FLM - PECULIARITA': strip da 12 tappi a profilo rientrante anti-interferenza vasta area ottica tipo extra-robusto
- FLM - IDONEE PER: biologia molecolare
- FLM - IDONEE PER: PCR
- FLM - IDONEE PER: PCR Real Time
- FLM - IDONEE PER: PCR Real Time FAST
- FLM - IDONEE PER: stoccaggio
- FLM - IDONEE PER: stoccaggio a -20°C

FLM - COMPATIBILE CON: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-SimpliAmp™-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
FLM - COMPATIBILE CON: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Well Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
FLM - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Deep Well
FLM - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Well Fast
FLM - COMPATIBILE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Thermomixer C-96 x 0.2ml
FLM - COMPATIBILE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-Opticon 2 System-96 x 0.2ml

FLM - COMPATIBILE CON: MJ Research Inc.-Opticon System-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
FLM - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Mx3000P™-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Mx3005P™-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Mx4000™-96 x 0.2ml
FLM - COMPATIBILE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TGradient-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TOptical-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-TRobot-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-UNO II-96 x 0.2ml
FLM - COMPATIBILE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml