

Microtubi da 0.2ml con tappo piatto ottico per PCR e PCR Real Time

Cf. 1000 pezzi - Microtubi da 0.2ml con tappo piatto ottico per PCR e PCR Real Time regular profile

Product codes:

Reference: 12 879 7601



Product gallery:



Product description:

Microtubi da 0,2ml, profilo regolare, eccellenti per PCR Real-Time

Compatibili con molti apparecchi per PCR Real-Time che accettano provettine da 0,2ml a profilo regolare.

Non compatibile con apparecchi FAST.

Cf. 1000pz.

COMPATIBILE CON

Applied Biosystems-310-96 x 0.2ml
Applied Biosystems-3100-96 x 0.2ml
Applied Biosystems-3130-96 x 0.2ml

Applied Biosystems-3700-96 x 0.2ml
Applied Biosystems-3730-96 x 0.2ml
Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
Applied Biosystems-ProFlex PCR System-32 x 0.2ml
Applied Biosystems-ProFlex PCR System-96 x 0.2ml
Applied Biosystems-QuantStudio 3 Real-Time PCR System-96 x 0.2ml
Applied Biosystems-QuantStudio 5 Real-Time PCR System-96 x 0.2ml
Applied Biosystems-QuantStudio 6 Flex-96 x 0.2ml
Applied Biosystems-QuantStudio 7 Flex-96 x 0.2ml
Applied Biosystems-QuantStudio 12k Flex-96 x 0.2ml
Applied Biosystems-SimpliAmp-96 x 0.2ml
Applied Biosystems-Veriti 96-Well Thermal Cycler-96 x 0.2ml
Applied Biosystems-Veriti Dx 96-Well Thermal Cycler-96 x 0.2ml
Applied Biosystems-ViiA 7 Real-Time PCR System-96 x 0.2ml
Bio-Rad Laboratories-C1000 Touch-96-Deep Well
Bio-Rad Laboratories-C1000 Touch-96-Well Fast
Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
Bio-Rad Laboratories-C1000-96-Deep Well
Bio-Rad Laboratories-C1000-96-Well Fast
Bio-Rad Laboratories-C1000-Dual 48/48 Fast
Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml

Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
Bio-Rad Laboratories-MyCycler-96 x 0.2ml
Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
Bio-Rad Laboratories-PTC-100-96 x 0.2ml
Bio-Rad Laboratories-S1000-96-Deep Well
Bio-Rad Laboratories-S1000-96-Well Fast
Bio-Rad Laboratories-S1000-Dual 48/48 Fast
Bio-Rad Laboratories-T100-96 x 0.2ml
Bio-Rad Laboratories-iCycler-48 x 0.2ml
Bio-Rad Laboratories-iCycler-96 x 0.2ml
Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
Eppendorf-Mastercycler-96 x 0.2ml
Eppendorf-Mastercycler-96-77 96 x 0.2ml
Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
Eppendorf-Mastercycler ep gradient-96 x 0.2ml
Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
Eppendorf-Mastercycler ep realplex-96 x 0.2ml
Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
Eppendorf-Mastercycler gradient-96 x 0.2ml
Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
Eppendorf-Mastercycler nexus-96 x 0.2ml
Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
Eppendorf-Mastercycler nexus eco-96 x 0.2ml
Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
Eppendorf-Mastercycler pro-96 x 0.2ml
Eppendorf-Mastercycler pro S-96 x 0.2ml
Eppendorf-Thermomixer C-96 x 0.2ml
Eppendorf-Thermomixer comfort-96 x 0.2ml
MJ Research Inc.-PTC-100-96 x 0.2ml
MJ Research Inc.-PTC-100-96 x 0.2ml AgV
MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml
MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml

MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
Roche-COBAS Taqman 48 Analyzer-24 x 0.2ml
Stratagene-Mx3000P-96 x 0.2ml
Stratagene-Mx3005P-96 x 0.2ml
Stratagene-Mx4000-96 x 0.2ml
Stratagene-Robocycler Gradient 96-96 x 0.2ml
Whatman Biometra-T1 Thermocycler-96 x 0.2ml
Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
Whatman Biometra-T3 Thermocycler-48 x 0.2ml
Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
Whatman Biometra-TGradient-96 x 0.2ml
Whatman Biometra-TOptical-96 x 0.2ml
Whatman Biometra-TPersonal-48 x 0.2ml
Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
Whatman Biometra-TProfessional-96 x 0.2ml
Whatman Biometra-TProfessional Basic-96 x 0.2ml
Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
Whatman Biometra-TProfessional TRIO-48 x 0.2ml
Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
Whatman Biometra-TRobot-96 x 0.2ml
Whatman Biometra-UNO-thermoblock-96 x 0.2ml
Whatman Biometra-UNO II-96 x 0.2ml

Product features:

MTB - TAPPO: piatto
MTB - CONFORMITA': non soggetto a CE
MTB - TIPO: PCR
MTB - VOLUME: 0.2ml
MTB - CAPACITA': 0.2ml
MTB - TAPPO FORABILE: no
MTB - TAPPO INSERITO: no
MTB - AREA OTTICA: 12.6mm
MTB - FORMA: regolare
MTB - FONDO: conico
MTB - SELF STANDING: no
MTB - MATERIALE: PP
MTB - GRADUATO: no
MTB - OPACITA': leggermente opaca
MTB - AREA SCRIVIBILE: no
MTB - RANGE TEMPERATURA: -30°C a +100°C
MTB - TEMPERATURA MAX: 100°C
MTB - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - AUTOCLAVABILE: si
MTB - STERILITA': no
MTB - DIMENSIONI: h 20.9/21.8mm
MTB - PECULIARITA': microtubi da 0.2ml con tappo piatto per PCR e PCR Real-Time
MTB - IDONEO PER: preparazione campioni
MTB - IDONEO PER: stoccaggio
MTB - IDONEO PER: PCR
MTB - IDONEO PER: PCR Real Time
MTB - RESISTENTE: congelamento a -20°C
MTB - PACKAGING: busta bulk
MTB - PARETI: sottili 0.3mm
MTB - EVAPORAZIONE CAMPIONE: minore del 3%
MTB - BINDING: DNA minore 0.3%
MTB - BINDING: proteine minore 0.8%
MTB - FORZA SEGNALE: 125
MTB - POSSIBILE INTERFERENZA SEGNALE: 70
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Deep Well
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Well Fast
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
MTB - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml

MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Deep Well
MTB - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Well Fast
MTB - COMPATIBILE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
MTB - COMPATIBILE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Thermomixer C-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml

MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TGradient-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TPersonal-48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TRobot-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-UNO II-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-310-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-3100-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-3130-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-3700-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-3730-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBILE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBILE CON: Applied Biosystems-ABI PRISM 7000 SDS-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-ABI PRISM 7700 SDS-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-ProFlex PCR System-32 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-ProFlex PCR System-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-QuantStudio 12k Flex-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-QuantStudio 3 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-QuantStudio 5 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-QuantStudio 6 Flex-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-QuantStudio 7 Flex-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-SimpliAmp-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-Veriti 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-Veriti Dx 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBILE CON: Applied Biosystems-ViiA 7 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml

MTB - COMPATIBILE CON: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
MTB - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBILE CON: Roche-COBAS Taqman 48 Analyzer-24 x 0.2ml
MTB - COMPATIBILE CON: Stratagene-Mx3000P™-96 x 0.2ml
MTB - COMPATIBILE CON: Stratagene-Mx3005P™-96 x 0.2ml
MTB - COMPATIBILE CON: Stratagene-Mx4000™-96 x 0.2ml
MTB - COMPATIBILE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
MTB - COMPATIBILE CON: Whatman Biometra-TOptical-96 x 0.2ml