

# Microtubi 0.2ml PCR tappo a cupola anti-contaminazione


Cf. 1000 pezzi - Microtubi da 0,2ml per PCR, tappo a cupola, anti-contaminazione



## Product codes:

Reference: 12 879 9511

## Product gallery:



## Product description:

Microtubi da 0,2ml, tappo a cupola, anti-contaminazione

Compatibili con molti apparecchi per PCR che accettano provettine da 0,2ml con tappo a cupola.

Profilo regolare.

Cf. 1000pz.

Visiti la scheda tecnica per consultare le caratteristiche eccellenti di questo prodotto!

COMPATIBILE CON:

MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	30-48 Alpha unit 48 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml
Whatman Biometra	T1 Thermocycler	96 x 0.2ml
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml
Whatman Biometra	T3 Thermocycler	48 x 0.2ml
Whatman Biometra	T3 Thermocycler Combi	48-18 48 x 0.2ml
Whatman Biometra	T3000 Thermocycler	48 x 0.2ml
Whatman Biometra	T3000 Thermocycler Combi	48-18 48 x 0.2ml
Whatman Biometra	TGradient	96 x 0.2ml
Whatman Biometra	TPersonal	48 x 0.2ml
Whatman Biometra	TPersonal Combi	48-18 48 x 0.2ml
Whatman Biometra	TProfessional	96 x 0.2ml
Whatman Biometra	TProfessional Basic	96 x 0.2ml
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml
Whatman Biometra	TProfessional Standard Gradient	96 x 0.2ml
Whatman Biometra	TProfessional TRIO	48 x 0.2ml
Whatman Biometra	TProfessional TRIO combi	48-30 30 x 0.5ml
Whatman Biometra	TProfessional TRIO combi	48-30 48 x 0.2ml
Whatman Biometra	TRobot	96 x 0.2ml
Whatman Biometra	UNO II	96 x 0.2ml
Whatman Biometra	UNO-thermoblock	96 x 0.2ml

### Product features:

MTB - TAPPO: semi-cupola  
 MTB - TAPPO: anticontaminazione  
 MTB - CONFORMITA': non soggetto a CE  
 MTB - TIPO: PCR  
 MTB - VOLUME: 0.2ml  
 MTB - CAPACITA': 0.2ml  
 MTB - TAPPO FORABILE: no  
 MTB - TAPPO INSERITO: no  
 MTB - AREA OTTICA: 6.8mm  
 MTB - FORMA: regolare  
 MTB - FONDO: conico  
 MTB - SELF STANDING: no  
 MTB - MATERIALE: PP  
 MTB - GRADUATO: no  
 MTB - OPACITA': ultra-trasparente  
 MTB - AREA SCRIVIBILE: no  
 MTB - RANGE TEMPERATURA: -30°C a +100°C  
 MTB - TEMPERATURA MAX: 100°C  
 MTB - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE  
 MTB - AUTOCLAVABILE: si  
 MTB - STERILITA': no

MTB - DIMENSIONI: h 20.9/21.8mm  
MTB - PECULIARITA': microtubi da 0.2ml con tappo a cupola anti-contaminazione  
MTB - IDONEO PER: preparazione campioni  
MTB - IDONEO PER: stoccaggio  
MTB - IDONEO PER: PCR  
MTB - IDONEO PER: PCR Real Time  
MTB - IDONEO PER: PCR Real Time FAST  
MTB - RESISTENTE: congelamento a -20°C  
MTB - PACKAGING: busta bulk  
MTB - PARETI: sottili 0.3mm  
MTB - EVAPORAZIONE CAMPIONE: minore del 3%  
MTB - BINDING: DNA minore 0.3%  
MTB - BINDING: proteine minore 0.8%  
MTB - FORZA SEGNALE: 100  
MTB - POSSIBILE INTERFERENZA SEGNALE: 100  
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml  
MTB - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Deep Well  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Well Fast  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Deep Well  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Well Fast  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast  
MTB - COMPATIBILE CON: Bio-Rad Laboratories-T100-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml

MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler pro-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Thermomixer C-96 x 0.2ml  
MTB - COMPATIBILE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml  
MTB - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TGradient-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TPersonal-48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml  
MTB - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-TRobot-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-UNO II-96 x 0.2ml  
MTB - COMPATIBILE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml