

# Piastra full-skirted 96x0.1ml (LP) Well Firm Deck Plate

Cf 20pz - Piastra full-skirted 96x0.1ml (LP) Well Firm Deck Plate - Ultra Trasparente



## Product codes:

Reference: 12 893 0701

## Product gallery:



## Product description:

Piastra da 0.1ml fully skirteds extra rigide, piatte, con bordo, impilabili e adatta alla robotica. possono essere utilizzate in diversi blocchi da 0,1 ml e 0,2 ml, termociclatori per PCR e Real-Time

Questa piastra FULL SKIRTED da 96 x 0,1 ml con piano rigido è progettata come alternativa diretta alle piastre a due componenti a 96 pozzetti con telaio rigido, come le piastre Bio-Rad Hard-Shell®, Eppendorf twin.tec® e Thermo FrameStar®

Le piastre Firm Deck sono estremamente piatte, rigide e indeformabili, progettate per resistere alle alte temperature durante la (q)PCR e la termosaldatura.

La chiusura può essere effettuata con qualsiasi striscia di tappi EU 8 o Opti-Seal™

Taglio dell'angolo a nord-est e sud-est

**COMPATIBILE CON**



|                                          |                                 |
|------------------------------------------|---------------------------------|
| Hain-LIFESCIENCE                         | GTQ-Cycler 96                   |
| Kyratec                                  | SuperCycler                     |
| LongGene Scientific Instruments Co.,Ltd. | A100                            |
| LongGene Scientific Instruments Co.,Ltd. | A200                            |
| LongGene Scientific Instruments Co.,Ltd. | A200                            |
| MJ Research Inc.                         | PTC-100                         |
| MJ Research Inc.                         | PTC-100                         |
| MJ Research Inc.                         | PTC-200 DNA Engine              |
| MJ Research Inc.                         | PTC-200 DNA Engine              |
| MJ Research Inc.                         | PTC-220 DNA Engine Dyad         |
| MJ Research Inc.                         | PTC-220 DNA Engine Dyad         |
| MJ Research Inc.                         | PTC-221 DNA Engine Disciple     |
| MJ Research Inc.                         | PTC-221 DNA Engine Disciple     |
| MJ Research Inc.                         | PTC-225 DNA Engine Tetrad       |
| MJ Research Inc.                         | PTC-225 DNA Engine Tetrad       |
| MJ Research Inc.                         | PTC-240 DNA Engine Tetrad 2     |
| MJ Research Inc.                         | PTC-240 DNA Engine Tetrad 2     |
| PCRMax                                   | AlphaCycler AC-1                |
| PCRMax                                   | AlphaCycler AC-2                |
| PCRMax                                   | AlphaCycler AC-4                |
| Qiagen                                   | QIAamplifier                    |
| Qiagen                                   | QIAquant 96 2plex               |
| Qiagen                                   | QIAquant 96 5plex               |
| Techne                                   | TC-4000                         |
| Techne                                   | TC-5000                         |
| Techne                                   | TC-Plus                         |
| Thermo Hybaid                            | MultiBlock System               |
| ThermoFisher Scientific                  | Arktik Thermal Cycler           |
| VWR                                      | DuoCycler                       |
| VWR                                      | DuoCycler fast                  |
| VWR                                      | QuattroCycler                   |
| VWR                                      | QuattroCycler fast              |
| VWR                                      | UnoCycler                       |
| VWR                                      | UnoCycler fast                  |
| Whatman Biometra                         | T1 Thermocycler                 |
| Whatman Biometra                         | TGradient                       |
| Whatman Biometra                         | TProfessional                   |
| Whatman Biometra                         | TProfessional Standard Gradient |
| Whatman Biometra                         | TRobot                          |

### Product features:

PSC - TIPO: piastre PCR  
PSC - VOLUME: 0.1ml  
PSC - FORMA: skirted  
PSC - APERTURA: normale  
PSC - POZZETTI: 96  
PSC - FONDO: a V  
PSC - CHIUSURA: tramite pellicola  
PSC - CHIUSURA: tramite tappi  
PSC - CHIUSURA: tappetino  
PSC - APPLICAZIONI: PCR

PSC - APPLICAZIONI: PCR Real-Time  
PSC - APPLICAZIONI: biologia molecolare  
PSC - AUTOCLAVABILE: sì  
PSC - COLORE: trasparente  
PSC - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE  
PSC - STERILITA': no  
PSC - MATERIALE: PP  
PSC - CONFORMITA': non soggetto a CE  
PSC - TEMPERATURA MAX: 100°C  
PSC - TEMPERATURA MINIMA: -30°C  
PSC - PACKAGE: bulk  
PSC - ACCESSORI OPZIONALI: strip di tappi  
PSC - ACCESSORI OPZIONALI: film sigillanti  
PSC - ACCESSORI OPZIONALI: tappetini sigillanti  
PSC - IDONEO PER: PCR  
PSC - IDONEO PER: PCR Real-Time  
PSC - IDONEO PER: stoccaggio  
PSC - IDONEO PER: preparazione campioni  
PSC - PROFILO: regolare  
PSC - COMPATIBILE CON: Agilent Technologies AriaDx Realtime PCR System 96 x 0.2ml  
PSC - COMPATIBILE CON: Agilent Technologies AriaMx Realtime PCR System 96 x 0.2ml  
PSC - COMPATIBILE CON: Analytik-Jena Alpha SC 96 x 0.2ml  
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler Monoblock 96 x 0.2ml  
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TOne 96 96 x 0.2ml  
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TOne 96G 96 x 0.2ml  
PSC - COMPATIBILE CON: Axygen® MaxyGene™ II 96 x 0.2ml  
PSC - COMPATIBILE CON: Bioer GenePro B-96G 96 x 0.2ml  
PSC - COMPATIBILE CON: Bioer LifePro TC-96/G/H(b)A 96 x 0.2ml  
PSC - COMPATIBILE CON: Bioer Line-Gene 9600 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Well Fast  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Well Fast  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus Deepwell 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC-100 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Well Fast  
PSC - COMPATIBILE CON: Blue-Ray Biotech TurboCycler TCST-9610 96 x 0.2ml  
PSC - COMPATIBILE CON: Blue-Ray Biotech TurboCycler TCST-9620 96 x 0.2ml  
PSC - COMPATIBILE CON: Bulldog Bio GenePro™ Thermal Cycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Bulldog Bio GeneTouch™ Thermal Cycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Bulldog Bio LifePro™ Thermal Cycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Bulldog Bio LifeTouch™ Thermal Cycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Corbett Research Palm-Cycler 96 x 0.2ml

PSC - COMPATIBILE CON: eENZYME GenePro Thermal Cyclers 96 x 0.2ml  
PSC - COMPATIBILE CON: eENZYME LifeExpress B-96D 96 x 0.2ml  
PSC - COMPATIBILE CON: eENZYME LifeExpress TC-96 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96-77 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro S 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50a 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50i 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50l 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50s 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50t 384 Microplate  
PSC - COMPATIBILE CON: Eppendorf Thermomixer C 96 x 0.2ml  
PSC - COMPATIBILE CON: Esco Swift™ Maxi 96 x 0.2ml  
PSC - COMPATIBILE CON: gradient  
PSC - COMPATIBILE CON: Esco Swift™ MaxPro 96 x 0.2ml  
PSC - COMPATIBILE CON: GE Healthcare MegaBACE 1000 96 x 0.1ml  
PSC - COMPATIBILE CON: GE Healthcare MegaBACE 4000 96 x 0.1ml  
PSC - COMPATIBILE CON: GE Healthcare MegaBACE 500 96 x 0.1ml  
PSC - COMPATIBILE CON: Hain-LIFESCIENCE GTQ-Cycler 96 96 x 0.2ml  
PSC - COMPATIBILE CON: Kyrtec SuperCycler 96 x 0.2ml  
PSC - COMPATIBILE CON: LongGene Scientific Instruments Co.  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml  
PSC - COMPATIBILE CON: AgV  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-1 96 x 0.2ml  
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-2 96 x 0.2ml  
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-4 96 x 0.2ml  
PSC - COMPATIBILE CON: Qiagen QIAamplifier 96 x 0.2ml  
PSC - COMPATIBILE CON: Qiagen QIAquant 96 2plex 96 x 0.2ml  
PSC - COMPATIBILE CON: Qiagen QIAquant 96 5plex 96 x 0.2ml  
PSC - COMPATIBILE CON: Techne TC-4000 96 x 0.2ml  
PSC - COMPATIBILE CON: Techne TC-5000 96 x 0.2ml  
PSC - COMPATIBILE CON: Techne TC-Plus 96 x 0.2ml  
PSC - COMPATIBILE CON: Thermo Hybrid MultiBlock System 96 x 0.2ml  
PSC - COMPATIBILE CON: ThermoFisher Scientific Arktik Thermal Cycler 96 x 0.2ml  
PSC - COMPATIBILE CON: VWR DuoCycler 96 x 0.2ml  
PSC - COMPATIBILE CON: VWR DuoCycler fast 96 x 0.2ml  
PSC - COMPATIBILE CON: VWR QuattroCycler 96 x 0.2ml  
PSC - COMPATIBILE CON: VWR QuattroCycler fast 96 x 0.2ml  
PSC - COMPATIBILE CON: VWR UnoCycler 96 x 0.2ml

PSC - COMPATIBILE CON: VWR UnoCycler fast 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TGradient 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TProfessional 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TRobot 96 x 0.2ml  
PSC - CAPACITA': 0.2ml  
PSC - PARETI: sottili 0.35mm  
PSC - TRATTAMENTO SEGNALE: pareti frosted  
PSC - FORZA SEGNALE: 125  
PSC - POSSIBILE INTERFERENZA SEGNALE: 70  
PSC - RESISTENTE: congelamento a -20°C  
PSC - EVAPORAZIONE CAMPIONE: minore del 3%  
PSC - BINDING: DNA minore 0.3%  
PSC - BINDING: proteine minore 0.8%  
PSC - GRIGLIA ALFANUMERICA: si