


# Piastra no-skirted 96x0.1ml (LP) Well Firm Deck Plate


Cf 40pz - Piastra Actares® No-skirted 96x0.1ml (LP) Well Firm Deck Plate- Ultra-trasp.



## Product codes:

Reference: 12 893 0761

## Product gallery:



## Product description:

Piastra da 0.1ml no skirted, 996 pozzetti, extra rigide, piatte, con bordo, impilabili sono progettate per applicazioni PCR e qPCR

Queste piastre sono compatibili con la maggior parte dei modelli di termociclatori (q)PCR che accettano prodotti a basso profilo (LP) rogettate per resistere alle alte temperature durante la (q)PCR e la termosaldatura.

La chiusura può essere effettuata con qualsiasi striscia di tappi EU 8 o Opti-Seal™

Taglio dell'angolo a nord-est e sud-est

**COMPATIBILE CON**



|                  |                                        |                                |
|------------------|----------------------------------------|--------------------------------|
| MJ Research Inc. | PTC-200 DNA Engine                     | Chromo4 96 x 0.2ml             |
| MJ Research Inc. | PTC-220 DNA Engine Dyad                | 96v Alpha unit 96 x 0.2ml      |
| MJ Research Inc. | PTC-220 DNA Engine Dyad                | 96v Moto Alpha Unit 96 x 0.2ml |
| MJ Research Inc. | PTC-220 DNA Engine Dyad                | Chromo4 96 x 0.2ml             |
| MJ Research Inc. | PTC-221 DNA Engine Disciple            | 96v Alpha unit 96 x 0.2ml      |
| MJ Research Inc. | PTC-221 DNA Engine Disciple            | 96v Moto Alpha Unit 96 x 0.2ml |
| MJ Research Inc. | PTC-221 DNA Engine Disciple            | Chromo4 96 x 0.2ml             |
| MJ Research Inc. | PTC-225 DNA Engine Tetrad              | 96v Alpha unit 96 x 0.2ml      |
| MJ Research Inc. | PTC-225 DNA Engine Tetrad              | 96v Moto Alpha Unit 96 x 0.2ml |
| MJ Research Inc. | PTC-225 DNA Engine Tetrad              | Chromo4 96 x 0.2ml             |
| MJ Research Inc. | PTC-240 DNA Engine Tetrad 2            | 96v Alpha unit 96 x 0.2ml      |
| MJ Research Inc. | PTC-240 DNA Engine Tetrad 2            | 96v Moto Alpha Unit 96 x 0.2ml |
| MJ Research Inc. | PTC-240 DNA Engine Tetrad 2            | Chromo4 96 x 0.2ml             |
| Roche            | LightCycler ® 480 Real-Time PCR System | 96 x 0.2ml                     |
| Roche            | LightCycler ® 96 Real-Time PCR System  | 96 x 0.2ml                     |
| Whatman Biometra | T1 Thermocycler                        | 96 x 0.2ml                     |
| Whatman Biometra | T1 Thermocycler Combi                  | 96-77 96 x 0.2ml               |
| Whatman Biometra | TGradient                              | 96 x 0.2ml                     |
| Whatman Biometra | TProfessional                          | 96 x 0.2ml                     |
| Whatman Biometra | TProfessional Basic                    | 96 x 0.2ml                     |
| Whatman Biometra | TProfessional Basic Gradient           | 96 x 0.2ml                     |
| Whatman Biometra | TProfessional Standard Gradient        | 96 x 0.2ml                     |
| Whatman Biometra | TRobot                                 | 96 x 0.2ml                     |
| Whatman Biometra | UNO II                                 | 96 x 0.2ml                     |
| Whatman Biometra | UNO-thermoblock                        | 96 x 0.2ml                     |

## Product features:

PSC - TIPO: piastre PCR  
 PSC - VOLUME: 0.1ml  
 PSC - FORMA: no-skirted  
 PSC - APERTURA: normale  
 PSC - POZZETTI: 96  
 PSC - FONDO: a V  
 PSC - CHIUSURA: tramite pellicola  
 PSC - CHIUSURA: tramite tappi  
 PSC - CHIUSURA: tappetino  
 PSC - APPLICAZIONI: PCR  
 PSC - APPLICAZIONI: PCR Real-Time  
 PSC - APPLICAZIONI: biologia molecolare  
 PSC - AUTOCLAVABILE: si  
 PSC - COLORE: trasparente  
 PSC - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE  
 PSC - STERILITA': no  
 PSC - MATERIALE: PP  
 PSC - CONFORMITA': non soggetto a CE  
 PSC - TEMPERATURA MAX: 100°C  
 PSC - TEMPERATURA MINIMA: -30°C  
 PSC - PACKAGE: bulk  
 PSC - ACCESSORI OPZIONALI: strip di tappi  
 PSC - ACCESSORI OPZIONALI: film sigillanti  
 PSC - ACCESSORI OPZIONALI: tappetini sigillanti

PSC - IDONEO PER: PCR  
PSC - IDONEO PER: PCR Real-Time  
PSC - IDONEO PER: stoccaggio  
PSC - IDONEO PER: preparazione campioni  
PSC - PROFILO: regolare  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Well Fast  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Well Fast  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC-100 96 x 0.2ml  
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Deep Well  
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Well Fast  
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96-77 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro S 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50a 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50i 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50l 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50s 96 x 0.2ml  
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50t 384 Microplate  
PSC - COMPATIBILE CON: Eppendorf Thermomixer C 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml  
PSC - COMPATIBILE CON: AgV  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml  
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TGradient 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TProfessional 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml  
PSC - COMPATIBILE CON: Whatman Biometra TRobot 96 x 0.2ml  
PSC - COMPATIBILE CON: Applied Biosystems 7500 Fast Real-Time PCR System 96 x 0.2ml

PSC - COMPATIBILE CON: Fast

PSC - COMPATIBILE CON: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml

PSC - COMPATIBILE CON: Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml

PSC - COMPATIBILE CON: Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml

PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems StepOnePlus™ Real-Time PCR System 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml

PSC - COMPATIBILE CON: Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml

PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX Duet Real-Time PCR System 96 x 0.2ml

PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: Bio-Rad Laboratories Opticon 2 96 x 0.2ml

PSC - COMPATIBILE CON: Eppendorf Mastercycler 96-77 77 x 0.5ml

PSC - COMPATIBILE CON: Eppendorf Thermomixer comfort 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. Opticon 2 System 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. Opticon System 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml

PSC - COMPATIBILE CON: Roche LightCycler® 480 Real-Time PCR System 96 x 0.2ml

PSC - COMPATIBILE CON: Roche LightCycler® 96 Real-Time PCR System 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra TProfessional Basic 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra UNO II 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra UNO-thermoblock 96 x 0.2ml

PSC - CAPACITA': 0.2ml

PSC - PARETI: sottili 0.35mm

PSC - TRATTAMENTO SEGNALE: pareti frosted

PSC - FORZA SEGNALE: 125

PSC - POSSIBILE INTERFERENZA SEGNALE: 70

PSC - RESISTENTE: congelamento a -20°C

PSC - EVAPORAZIONE CAMPIONE: minore del 3%

PSC - BINDING: DNA minore 0.3%

PSC - BINDING: proteine minore 0.8%

PSC - GRIGLIA ALFANUMERICA: si