



Piastra Tear-Off da 12 strip da 8x0.1ml separabili

Cf 36pz - Piastra Tear-Off Actares® fully-skirted composta da 12 strip da 8x0.1ml (LP) Well Firm Deck separabili- Ultra-trasp

**Product codes:**

Reference: 12 893 0771

Product gallery:**Product description:**

Piastra Tear-Off Actares® fully-skirted composta da 12 strip da 8x0.1ml (LP) Well Firm Deck separabili- Ultra-trasp

Compatibile con quasi tutti i modelli di termociclatori PCR e qPCR (fast) che accettano prodotti a basso profilo come termociclatori ABI/LT/Thermo Fast, Roche®, Bio-Rad®, Eppendorf® e altri.

Progettato per applicazioni PCR e qPCR e piattaforme di manipolazione dei liquidi

COMPATIBILE CON:

TaKaRa Bio Inc.	PCR Thermal Cycler Dice TP990	96 x 0.1ml
ThermoFisher Scientific	Arktik Thermal Cycler	48 x 0.2ml
VWR	DuoCycler	96 x 0.2ml
VWR	DuoCycler fast	96 x 0.2ml
VWR	QuattroCycler	96 x 0.2ml
VWR	QuattroCycler fast	96 x 0.2ml
VWR	UnoCycler	96 x 0.2ml
VWR	UnoCycler fast	96 x 0.2ml
Whatman Biometra	T1 Thermocycler	96 x 0.2ml
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml
Whatman Biometra	T3 Thermocycler	48 x 0.2ml
Whatman Biometra	T3 Thermocycler Combi	48-18 48 x 0.2ml
Whatman Biometra	T3000 Thermocycler	48 x 0.2ml
Whatman Biometra	T3000 Thermocycler Combi	48-18 48 x 0.2ml
Whatman Biometra	TGradient	96 x 0.2ml
Whatman Biometra	TOptical	96 x 0.2ml
Whatman Biometra	TPersonal	48 x 0.2ml
Whatman Biometra	TPersonal Combi	48-18 48 x 0.2ml
Whatman Biometra	TProfessional	96 x 0.2ml
Whatman Biometra	TProfessional Basic	96 x 0.2ml
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml
Whatman Biometra	TProfessional Standard Gradient	96 x 0.2ml
Whatman Biometra	TProfessional TRIO	48 x 0.2ml
Whatman Biometra	TProfessional TRIO combi	48-30 48 x 0.2ml
Whatman Biometra	TRobot	96 x 0.2ml
Whatman Biometra	UNO II	96 x 0.2ml
Whatman Biometra	UNO-thermoblock	96 x 0.2ml
Xian Tianlong Science & Technology	DTC-4	96-77 96 x 0.2ml
Xian Tianlong Science & Technology	DTC-4G	96-77 96 x 0.2ml

Product features:

PSC - TIPO: piastre PCR
 PSC - VOLUME: 0.1ml
 PSC - FORMA: semi-skirted
 PSC - APERTURA: normale
 PSC - POZZETTI: 96
 PSC - FONDO: a V
 PSC - CHIUSURA: tramite pellicola
 PSC - CHIUSURA: tramite tappi
 PSC - CHIUSURA: tappetino
 PSC - APPLICAZIONI: PCR
 PSC - APPLICAZIONI: PCR Real-Time FAST
 PSC - APPLICAZIONI: biologia molecolare
 PSC - AUTOCLAVABILE: sì
 PSC - COLORE: trasparente
 PSC - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE
 PSC - STERILITA': no
 PSC - MATERIALE: PP
 PSC - CONFORMITA': non soggetto a CE

PSC - TEMPERATURA MAX: 100°C
PSC - TEMPERATURA MINIMA: -30°C
PSC - PACKAGE: bulk
PSC - ACCESSORI OPZIONALI: strip di tappi
PSC - ACCESSORI OPZIONALI: film sigillanti
PSC - ACCESSORI OPZIONALI: tappetini sigillanti
PSC - IDONEO PER: PCR
PSC - IDONEO PER: PCR Real-Time
PSC - IDONEO PER: stoccaggio
PSC - IDONEO PER: preparazione campioni
PSC - PROFILO: regolare
PSC - COMPATIBILE CON: Agilent Technologies AriaDx Realtime PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Agilent Technologies AriaMx Realtime PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena Alpha SC 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler Monoblock 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TOne 96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TOne 96G 96 x 0.2ml
PSC - COMPATIBILE CON: Axygen® MaxyGene™ II 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer LifePro TC-96/G/H(b)A 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus Deepwell 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBILE CON: Blue-Ray Biotech TurboCycler TCST-9610 96 x 0.2ml
PSC - COMPATIBILE CON: Blue-Ray Biotech TurboCycler TCST-9620 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GenePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GeneTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio LifePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio LifeTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME GenePro Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME LifeExpress B-96D 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME LifeExpress TC-96 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50a 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50i 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50l 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50s 96 x 0.2ml

PSC - COMPATIBILE CON: Eppendorf Mastercycler X50t 384 Microplate
PSC - COMPATIBILE CON: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ Maxi 96 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ MaxPro 96 x 0.2ml
PSC - COMPATIBILE CON: Hain-LIFESCIENCE GTQ-Cycler 96 96 x 0.2ml
PSC - COMPATIBILE CON: Kyrtec SuperCycler 96 x 0.2ml
PSC - COMPATIBILE CON: LongGene Scientific Instruments Co.
PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-1 96 x 0.2ml
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-2 96 x 0.2ml
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-4 96 x 0.2ml
PSC - COMPATIBILE CON: Qiagen QIAamplifier 96 x 0.2ml
PSC - COMPATIBILE CON: Qiagen QIAquant 96 2plex 96 x 0.2ml
PSC - COMPATIBILE CON: Qiagen QIAquant 96 5plex 96 x 0.2ml
PSC - COMPATIBILE CON: VWR DuoCycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR DuoCycler fast 96 x 0.2ml
PSC - COMPATIBILE CON: VWR QuattroCycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR QuattroCycler fast 96 x 0.2ml
PSC - COMPATIBILE CON: VWR UnoCycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR UnoCycler fast 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems StepOnePlus™ Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX Duet Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories Opticon 2 96 x 0.2ml
PSC - COMPATIBILE CON: Roche LightCycler ® 480 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Roche LightCycler ® 96 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra UNO II 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra UNO-thermoblock 96 x 0.2ml
PSC - COMPATIBILE CON: Acugen Systems AG-9600 thermal station 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena Alpha SC 48 x 0.2ml

PSC - COMPATIBILE CON: Analytik-Jena Alpha SC 48-96 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler Twinblock 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler Twinblock Mix 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler2 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler2 48-18 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler2 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 G 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 SG 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced Twin 48 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced Twin 48 G 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic Gradient XL96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic XL96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard Gradient SL 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard SL 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional TRIO 48 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional TRIO Combi 48-18 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TRIO 48 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TRIO Combi 48-18 48 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TRobot II 96 G 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TRobot II 96 SG 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC201 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC401 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC701 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC801 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 6 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 7 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems StepOne™ Real-Time PCR System 48 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems ViiA™ 7 Real-Time PCR System Fast 96 x 0.1ml
PSC - COMPATIBILE CON: Barnstead Thermolyne Amplitron I 96 x 0.2ml
PSC - COMPATIBILE CON: Barnstead Thermolyne Amplitron II 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneExplorer GE-48D 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneExplorer GE-96G 96 x
PSC - COMPATIBILE CON: 0.2ml
PSC - COMPATIBILE CON: Bioer GenePro B-3048 30-48 30 x 0.5ml
PSC - COMPATIBILE CON: Bioer GenePro B-3048 30-48 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer GenePro B-48D 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer
PSC - COMPATIBILE CON: GenePro B-96G 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneQ TC-24/H(b) 24 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneTouch™ Thermal Cycler 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneTouch™ Thermal Cycler 96 x 0.2ml

PSC - COMPATIBILE CON: Bioer Life ECO Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer Life Express TC-48/H(t) 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer Life Express TC-96/T/H(a) 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer LifeTouch TC-96/G/H(b)B 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer Little Genius TC-24/H 24 x 0.2ml
PSC - COMPATIBILE CON: Bioer XPcycler TC-XP-A 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer XPcycler TC-XP-D 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer XPcycler TC-XP-E 48-30 48 x 0.2ml
PSC - COMPATIBILE CON: Bioer XPcycler TC-XP-G 96 x 0.2ml
PSC - COMPATIBILE CON: Biomérieux GENE-UP 96 x 0.2ml Fast
PSC - COMPATIBILE CON: Bioneer Exicycler 96 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Dual 48/48 Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch Dual 48/48 Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 6-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories MJ Mini 48-12 48 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC Tempo 96 Thermal Cyclers 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 Dual 48/48 Fast
PSC - COMPATIBILE CON: Biosan Ltd. TS-100 Thermo-Shaker 96 x 0.2ml
PSC - COMPATIBILE CON: Boeco TC-PRO 3048 30-48 48 x 0.2ml
PSC - COMPATIBILE CON: Boeco TC-PRO 48D 48 x 0.2ml
PSC - COMPATIBILE CON: Boeco TC-PRO 96G 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GenePro™ Thermal Cyclers 30-48 48 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GenePro™ Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GeneQ™ Thermal Cyclers 24 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GeneTouch™ Thermal Cyclers 30-48 48 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GeneTouch™ Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBILE CON: Chai Biotechnologies Inc. OpenQPCR Dual Channel 16 x 0.2ml
PSC - COMPATIBILE CON: hai Biotechnologies Inc. OpenQPCR Single Channel 8 x 0.2ml
PSC - COMPATIBILE CON: DuPont Qualicon BAX System Q7 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME GenePro Thermal Cyclers 48 x 0.2ml
PSC - COMPATIBILE CON: eENZYME GeneQ Personal Thermal Cyclers 24 x 0.2ml
PSC - COMPATIBILE CON: eENZYME LifeExpress TC-48 48 x 0.2ml
PSC - COMPATIBILE CON: eENZYME MyGene MG-96 Gradient Thermal Cyclers 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME MyGene MG-96 Thermal Cyclers 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler gradient 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus eco 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml
PSC - COMPATIBILE CON: Ericom Deltacycler I 96 x 0.2ml
PSC - COMPATIBILE CON: Ericom Deltacycler II 96 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ Maxi 48 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ Maxi 48-30 48 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ Maxi 96 x 0.2ml gradient
PSC - COMPATIBILE CON: Esco Swift™ MaxPro 48-30 48 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ MaxPro 48-48 48 x 0.2ml
PSC - COMPATIBILE CON: Esco Swift™ MiniPro 24 x 0.2ml
PSC - COMPATIBILE CON: EuroClone One Advanced 96 x 0.2ml

PSC - COMPATIBILE CON: EuroClone One Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: EuroClone peqSTAR 96 HPL 96 x 0.2ml
PSC - COMPATIBILE CON: EuroClone peqSTAR 96 ML 96 x 0.2ml
PSC - COMPATIBILE CON: EuroClone peqSTAR 96 Universal 96 x 0.2ml
PSC - COMPATIBILE CON: Finnzymes Instruments Piko™ 24-well 24-wells
PSC - COMPATIBILE CON: Finnzymes Instruments TCA0096 Arktik 96 x 0.2ml
PSC - COMPATIBILE CON: Finnzymes Instruments TCA4848 Arktik 2 x 48 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS1 48 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS1 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS1 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS2 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS4 48 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS4 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS4 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS4822 48 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GSX 48 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GSX 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GSX 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Hangzhou Jingle Scientific Instrument Co.
PSC - COMPATIBILE CON: Helena Biosciences Phoenix Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Lab-Line Instruments Programmable Thermal Block 30-48 48 x 0.2ml
PSC - COMPATIBILE CON: Lab-Line Instruments Programmable Thermal Block 96 x 0.2ml
PSC - COMPATIBILE CON: Labnet MultiGene™ Mini Personal Thermal Cycler 24 x 0.2ml
PSC - COMPATIBILE CON: Labnet MultiGene™ Thermal Cycler Gradient TC9600-G 96 x 0.2ml
PSC - COMPATIBILE CON: Labnet MultiGene™ Thermal Cycler TC9600 96 x 0.2ml
PSC - COMPATIBILE CON: Landgraf Varius M 96 x 0.2ml
PSC - COMPATIBILE CON: Landgraf Varius V 80 16 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml AgV
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 30 x 0.5ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml
PSC - COMPATIBILE CON: SensoQuest Labcycler Thermoblock 96
PSC - COMPATIBILE CON: TaKaRa Bio Inc. PCR Thermal Cycler Dice TP950 96 x 0.1ml
PSC - COMPATIBILE CON: TaKaRa Bio Inc. PCR Thermal Cycler Dice TP990 96 x 0.1ml
PSC - COMPATIBILE CON: ThermoFisher Scientific Arktik Thermal Cycler 48 x 2ml
PSC - COMPATIBILE CON: Whatman Biometra T3 Thermocycler 48 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T3 Thermocycler Combi 48-18 48 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T3000 Thermocycler 48 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T3000 Thermocycler Combi 48-18 48 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TOptical 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TPersonal 48 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TPersonal Combi 48-18 48 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra TProfessional TRIO 48 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional TRIO combi 48-30 48 x 0.2ml
PSC - COMPATIBILE CON: Xian Tianlong Science & Technology DTC-4 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Xian Tianlong Science & Technology DTC-4G 96-77 96 x 0.2ml
PSC - CAPACITA': 0.1ml
PSC - PARETI: sottili 0.35mm
PSC - TRATTAMENTO SEGNALE: pareti frosted
PSC - FORZA SEGNALE: 125
PSC - POSSIBILE INTERFERENZA SEGNALE: 70
PSC - RESISTENTE: congelamento a -20°C
PSC - EVAPORAZIONE CAMPIONE: minore del 3%
PSC - BINDING: DNA minore 0.3%
PSC - BINDING: proteine minore 0.8%
PSC - GRIGLIA ALFANUMERICA: si