



Piastra semi-skirted 96x0.1ml - trasparente

Cf 36pz - Piastra Actares® semi-skirted 96x0.1ml (LP) Well Firm Deck Plate- Design Universale - Ultra-trasp

**Product codes:**

Reference: 12 893 0801

Product gallery:**Product description:**

Piastra Actares® semi-skirted 96x0.1ml (LP) Well Firm Deck Plate- Design Universale - Ultra-trasp.

Piastra FD2 (Firm Deck) robusta da 96 x 0,1 ml, adatta a tutti i termociclatori qPCR, semi-skirted (DR8) Quindi non più 3 piastre diverse con angoli tagliati specifici per termociclatori in tempo reale Thermo, Roche, Bio-Rad e Qiagen, ma 1 piastra con 3 angoli tagliati per tutti i termociclatori.

progettate come alternativa diretta alle piastre a due componenti a 96 pozzetti con struttura rigida. Le piastre FD2 sono estremamente piatte, rigide e indeformabili, progettate per resistere alle alte temperature durante la qPCR e la termosaldatura

Taglio degli angoli a nord-ovest, nord-est e sud-est

COMPATIBILE CON:

Techne	Genius	96 x 0.2ml
Techne	Prime thermal cycler	96 x 0.2ml
Techne	TC-4000	96 x 0.2ml
Techne	TC-412	96 x 0.2ml
Techne	TC-5000	96 x 0.2ml
Techne	TC-512	96 x 0.2ml
Techne	TC-Plus	96 x 0.2ml
Techne	Touchgene Gradient	96 x 0.2ml
ThermoFisher Scientific	Arktik Thermal Cycler	96 x 0.2ml
VWR	DuoCycler	96 x 0.2ml
VWR	DuoCycler fast	96 x 0.2ml
VWR	QuattroCycler	96 x 0.2ml
VWR	QuattroCycler fast	96 x 0.2ml
VWR	UnoCycler	96 x 0.2ml
VWR	UnoCycler fast	96 x 0.2ml
Whatman Biometra	T1 Thermocycler	96 x 0.2ml
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml
Whatman Biometra	TGradient	96 x 0.2ml
Whatman Biometra	TProfessional	96 x 0.2ml
Whatman Biometra	TProfessional Basic	96 x 0.2ml
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml
Whatman Biometra	TProfessional Standard Gradient	96 x 0.2ml
Whatman Biometra	TRobot	96 x 0.2ml
Whatman Biometra	UNO II	96 x 0.2ml
Whatman Biometra	UNO-thermoblock	96 x 0.2ml

Product features:

PSC - TIPO: piastre PCR
 PSC - VOLUME: 0.1ml
 PSC - FORMA: semi-skirted
 PSC - APERTURA: normale
 PSC - POZZETTI: 96
 PSC - FONDO: a V
 PSC - CHIUSURA: tramite pellicola
 PSC - CHIUSURA: tramite tappi
 PSC - CHIUSURA: tappetino
 PSC - APPLICAZIONI: PCR
 PSC - APPLICAZIONI: PCR Real-Time FAST
 PSC - APPLICAZIONI: biologia molecolare
 PSC - AUTOCLAVABILE: si
 PSC - COLORE: trasparente
 PSC - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE
 PSC - STERILITA': no
 PSC - MATERIALE: PP
 PSC - CONFORMITA': non soggetto a CE
 PSC - TEMPERATURA MAX: 100°C
 PSC - TEMPERATURA MINIMA: -30°C
 PSC - PACKAGE: bulk
 PSC - ACCESSORI OPZIONALI: strip di tappi
 PSC - ACCESSORI OPZIONALI: film sigillanti
 PSC - ACCESSORI OPZIONALI: tappetini sigillanti

PSC - IDONEO PER: PCR
PSC - IDONEO PER: PCR Real-Time
PSC - IDONEO PER: stoccaggio
PSC - IDONEO PER: preparazione campioni
PSC - PROFILO: regolare
PSC - COMPATIBILE CON: Analytik-Jena Alpha SC 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler Monoblock 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TOne 96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TOne 96G 96 x 0.2ml
PSC - COMPATIBILE CON: Axygen® MaxyGene™ II 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer GenePro B-96G 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer LifePro TC-96/G/H(b)A 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer Line-Gene 9600 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Opus Deepwell 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBILE CON: Blue-Ray Biotech TurboCycler TCST-9610 96 x 0.2ml
PSC - COMPATIBILE CON: Blue-Ray Biotech TurboCycler TCST-9620 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GenePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio GeneTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio LifePro™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bulldog Bio LifeTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME GenePro Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME LifeExpress B-96D 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME LifeExpress TC-96 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50a 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50i 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50l 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50s 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler X50t 384 Microplate
PSC - COMPATIBILE CON: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBILE CON: Hain-LIFESCIENCE GTQ-Cycler 96 96 x 0.2ml
PSC - COMPATIBILE CON: Kyrtec SuperCycler 96 x 0.2ml
PSC - COMPATIBILE CON: LongGene Scientific Instruments Co.

PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-1 96 x 0.2ml
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-2 96 x 0.2ml
PSC - COMPATIBILE CON: PCRMax AlphaCycler AC-4 96 x 0.2ml
PSC - COMPATIBILE CON: Qiagen QIAamplifier 96 x 0.2ml
PSC - COMPATIBILE CON: Qiagen QIAquant 96 2plex 96 x 0.2ml
PSC - COMPATIBILE CON: Qiagen QIAquant 96 5plex 96 x 0.2ml
PSC - COMPATIBILE CON: Techne TC-4000 96 x 0.2ml
PSC - COMPATIBILE CON: Techne TC-5000 96 x 0.2ml
PSC - COMPATIBILE CON: Techne TC-Plus 96 x 0.2ml
PSC - COMPATIBILE CON: ThermoFisher Scientific Arktik Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR DuoCycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR DuoCycler fast 96 x 0.2ml
PSC - COMPATIBILE CON: VWR QuattroCycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR QuattroCycler fast 96 x 0.2ml
PSC - COMPATIBILE CON: VWR UnoCycler 96 x 0.2ml
PSC - COMPATIBILE CON: VWR UnoCycler fast 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TGradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems StepOnePlus™ Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX Duet Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories Opticon 2 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. Opticon 2 System 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. Opticon System 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Roche LightCycler® 480 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Roche LightCycler® 96 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra UNO II 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra UNO-thermoblock 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena FlexCycler2 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 G 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 SG 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic Gradient XL96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Basic XL96 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard Gradient SL 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TProfessional Standard SL 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC201 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC401 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC701 96 x 0.2ml
PSC - COMPATIBILE CON: Apollo ATC801 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 6 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml Fast
PSC - COMPATIBILE CON: Applied Biosystems QuantStudio™ 7 PRO Real-Time PCR System 96 x 0.1ml
PSC - COMPATIBILE CON: Applied Biosystems ViiA™ 7 Real-Time PCR System Fast 96 x 0.1ml
PSC - COMPATIBILE CON: Barnstead Thermolyne Amplitron I 96 x 0.2ml
PSC - COMPATIBILE CON: Barnstead Thermolyne Amplitron II 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneTouch™ Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer Life Express TC-96/T/H(a) 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer LifeTouch TC-96/G/H(b)B 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer XPcycler TC-XP-A 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC Tempo 96 Thermal Cycler 96 x 0.2ml
PSC - COMPATIBILE CON: eENZYME MyGene MG-96 Gradient Thermal Cycler 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler gradient 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus eco 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml
PSC - COMPATIBILE CON: Finnzymes Instruments TCA0096 Arktik 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS1 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS1 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS2 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS4 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GS4 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GSX 96 x 0.2ml
PSC - COMPATIBILE CON: Gene Technologies Ltd. G-Storm GSX 96-48 96 x 0.2ml
PSC - COMPATIBILE CON: SensoQuest Labcycler Thermoblock 96
PSC - COMPATIBILE CON: TaKaRa Bio Inc. PCR Thermal Cycler Dice TP950 96 x 0.1ml
PSC - COMPATIBILE CON: TaKaRa Bio Inc. PCR Thermal Cycler Dice TP990 96 x 0.1ml
PSC - COMPATIBILE CON: Analytik-Jena qTower 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena qTOWER³ touch 96 x 0.2ml

PSC - COMPATIBILE CON: Analytik-Jena SpeedCycler² 96 x 0.2ml
PSC - COMPATIBILE CON: Analytik-Jena - Biometra TAdvanced 96 S 96 x 0.2ml
PSC - COMPATIBILE CON: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Standard
PSC - COMPATIBILE CON: Azure biosystems Cielo 3 Real Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Azure biosystems Cielo 6 Real Time PCR System 96 x 0.2ml
PSC - COMPATIBILE CON: Bioer GeneExplorer GE-96G 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC Tempo Deepwell Thermal Cycler 96-Deep Well Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBILE CON: Biosystems MJ96+ 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler gradient 96-77 77 x 0.5ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus GSX1 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus GSX1e 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus SX1 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler nexus SX1e 96 x 0.2ml
PSC - COMPATIBILE CON: Landgraf Varius V 45 9 x 0.5ml
PSC - COMPATIBILE CON: AgV MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Peqlab Primus 96 Advanced 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Roche cobas z 480 96 x 0.2ml
PSC - COMPATIBILE CON: Seegene Inc. SEEAMP™ 96 x 0.2ml
PSC - COMPATIBILE CON: Techne Flexigene 96 x 0.2ml
PSC - COMPATIBILE CON: Techne Gene-E 96 x 0.2ml
PSC - COMPATIBILE CON: Techne Genius 96 x 0.2ml
PSC - COMPATIBILE CON: Techne Prime thermal cycler 96 x 0.2ml
PSC - COMPATIBILE CON: Techne TC-412 96 x 0.2ml
PSC - COMPATIBILE CON: Techne TC-512 96 x 0.2ml
PSC - COMPATIBILE CON: Techne Touchgene Gradient 96 x 0.2ml
PSC - CAPACITA': 0.1ml
PSC - PARETI: sottili 0.35mm
PSC - TRATTAMENTO SEGNALE: pareti frosted
PSC - FORZA SEGNALE: 125
PSC - POSSIBILE INTERFERENZA SEGNALE: 70
PSC - RESISTENTE: congelamento a -20°C
PSC - EVAPORAZIONE CAMPIONE: minore del 3%
PSC - BINDING: DNA minore 0.3%
PSC - BINDING: proteine minore 0.8%
PSC - GRIGLIA ALFANUMERICA: si