

Strip 8x0.1ml extra robuste bianche

strip 8x0,1ml ideale per apparecchi tipo FAST o a basso profilo,

Extra robuste - bianche

Cf. 120 strip

Product codes:

Reference: 12 954 8391



Product gallery:



Product description:

Strips da 0,1ml, profilo basso, eccellenti per PCR Real-Time.

Tutta la serie di strips possiede una forma asimmetrica come riferimento orientativo, un'area scrivibile e dei fori per fissarli sulla griglia atta ad ospitare le strips per creare la propria piastra.

Extra-robuste e bianche (+ 90% segnale), pareti estremamente uniformi e sottili.

Prodotte con una speciale miscela in PP studiata per applicazioni in PCR e qPCR,

basso grado di evaporazione (inferiore al 3%, idoneo anche per volumi minori di 5µl)

Low binding di proteine

DNA, DNase, RNA, RNase & Pyrogen FREE

Autocavabili

Cf. 120pz.

COMPATIBILI CON:

Applied Biosystems 7500 Fast Real-Time PCR System 96 x 0.2ml Fast Applied Biosystems 9800 Fast Thermal Cycler 96 x 0.2ml Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml Applied Biosystems QuantStudio™ 12k Flex 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 6 Flex 96 x 0.1ml Fast Applied Biosystems QuantStudio™ 7 Flex 96 x 0.1ml Fast Applied Biosystems StepOne™ Real-Time PCR System 48 x 0.1ml Applied Biosystems Veriti™ 96-Well Fast Thermal Cycler 96 x 0.1ml Applied Biosystems Veriti™ Dx 96-Well Fast Thermal Cycler 96 x 0.1ml Applied Biosystems ViiA™ 7 Real-Time PCR System Fast 96 x 0.1ml Bio-Rad Laboratories C1000 96-Well Fast Bio-Rad Laboratories C1000 Dual 48/48 Fast Bio-Rad Laboratories C1000 Touch 96-Deep Well Bio-Rad Laboratories C1000 Touch 96-Well Fast Bio-Rad Laboratories C1000 Touch Dual 48/48 Fast Bio-Rad Laboratories DNA Engine 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Disciple 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml Bio-Rad Laboratories DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml Bio-Rad Laboratories MJ Mini 48-12 48 x 0.2ml Bio-Rad Laboratories PTC-100 96 x 0.2ml Bio-Rad Laboratories S1000 96-Deep Well Bio-Rad Laboratories S1000 96-Well Fast Bio-Rad Laboratories S1000 Dual 48/48 Fast Eppendorf Mastercycler 96 x 0.2ml Eppendorf Mastercycler 96-77 96 x 0.2ml Eppendorf Mastercycler ep gradient 96 x 0.2ml Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml Eppendorf Mastercycler gradient 96 x 0.2ml Eppendorf Mastercycler gradient 96-77 96 x 0.2ml Eppendorf Mastercycler nexus 96 x 0.2ml Eppendorf Mastercycler nexus eco 96 x 0.2ml Eppendorf Mastercycler nexus gradient 96 x 0.2ml Eppendorf Mastercycler nexus gradient eco 96 x 0.2ml Eppendorf Mastercycler pro 96 x 0.2ml Eppendorf Mastercycler pro S 96 x 0.2ml Eppendorf Thermomixer C 96 x 0.2ml MJ Research Inc. PTC-100 96 x 0.2ml MJ Research Inc. PTC-100 96 x 0.2ml AgV MJ Research Inc. PTC-200 DNA Engine 30-48 Alpha unit 30 x 0.5ml MJ Research Inc. PTC-200 DNA Engine 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-225

DNA Engine Tetrad 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 30-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 48-48 Alpha unit 48 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml Roche LightCycler ® Nano 32 x 0.2ml Whatman Biometra T1 Thermocycler 96 x 0.2ml Whatman Biometra T1 Thermocycler Combi 96-77 96 x 0.2ml Whatman Biometra T3 Thermocycler 48 x 0.2ml Whatman Biometra T3 Thermocycler Combi 48-18 48 x 0.2ml Whatman Biometra T3000 Thermocycler 48 x 0.2ml Whatman Biometra T3000 Thermocycler Combi 48-18 48 x 0.2ml Whatman Biometra TGradient 96 x 0.2ml Whatman Biometra TOptical 96 x 0.2ml Whatman Biometra TPersonal 48 x 0.2ml Whatman Biometra TPersonal Combi 48-18 48 x 0.2ml Whatman Biometra TProfessional 96 x 0.2ml Whatman Biometra TProfessional Basic 96 x 0.2ml Whatman Biometra TProfessional Basic Gradient 96 x 0.2ml Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml Whatman Biometra TProfessional TRIO 48 x 0.2ml Whatman Biometra TProfessional TRIO combi 48-30 48 x 0.2ml Whatman Biometra TRobot 96 x 0.2ml Whatman Biometra UNO II 96 x 0.2ml Whatman Biometra UNO-thermoblock 96 x 0.2ml

Product features:

SP - TIPO: strip da 8 microtubi
SP - TAPPO: senza
SP - VOLUME: 0.1ml
SP - CAPACITA': 0.2ml
SP - POZZETTI: 8
SP - PROFILO: basso tipo FAST/0.1
SP - ALTEZZA: 15.6 / 16.6mm ca
SP - FORMA: striscia robusta
SP - FONDO: conico
SP - MATERIALE: PP
SP - APPLICAZIONI: PCR
SP - APPLICAZIONI: PCR Real Time
SP - APPLICAZIONI: PCR Real Time FAST
SP - AUTOCLAVABILE: si
SP - COLORE: bianco
SP - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITA': no
SP - CONFORMITA': non soggetto a CE
SP - TEMPERATURA MAX: 100°
SP - TEMPERATURA MINIMA: -30°C
SP - PACKAGE: busta
SP - PECULIARITA': strip da 8x0.1 per PCR Real Time FAST e PCR extra robuste bianco
SP - ACCESSORI OPZIONALI: strip di tappi
SP - ACCESSORI OPZIONALI: griglia assemblaggio bianca
SP - IDONEO PER: PCR
SP - IDONEO PER: PCR Real Time

SP - IDONEO PER: PCR Real Time FAST
SP - IDONEO PER: stoccaggio
SP - PARETI: sottili 0.30mm
SP - TRATTAMENTO SEGNALE: pareti bianche
SP - FORZA SEGNALE: 180
SP - POSSIBILE INTERFERENZA SEGNALE: inferiore a 2
SP - RESISTENTE: congelamento a -20°C
SP - EVAPORAZIONE CAMPIONE: minore del 3%
SP - BINDING: DNA minore 0.3%
SP - BINDING: proteine minore 0.8%
SP - COMPATIBILE CON: Applied Biosystems-7500 Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBILE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Fast
SP - COMPATIBILE CON: Applied Biosystems-9800 Fast Thermal Cycler-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-QuantStudio 12k Flex-96 x 0.1ml Fast
SP - COMPATIBILE CON: Applied Biosystems-QuantStudio 3 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBILE CON: Applied Biosystems-QuantStudio 5 Real-Time PCR System-96 x 0.1ml Fast
SP - COMPATIBILE CON: Applied Biosystems-QuantStudio 6 Flex-96 x 0.1ml Fast
SP - COMPATIBILE CON: Applied Biosystems-QuantStudio 7 Flex-96 x 0.1ml Fast
SP - COMPATIBILE CON: Applied Biosystems-StepOnePlus Real-Time PCR System-96 x 0.1ml
SP - COMPATIBILE CON: Applied Biosystems-Veriti 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBILE CON: Applied Biosystems-Veriti Dx 96-Well Fast Thermal Cycler-96 x 0.1ml
SP - COMPATIBILE CON: Applied Biosystems-ViiA 7 Real-Time PCR System Fast-96 x 0.1ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-Opticon 2-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-Opticon 2 System-96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-Opticon System-96 x 0.2ml

SP - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
SP - COMPATIBILE CON: Roche-LightCycler ® 96 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-StepOne Real-Time PCR System-48 x 0.1ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Roche-LightCycler ® 480 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml

SP - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml

SP - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml

SP - COMPATIBILE CON: Whatman Biometra-TPersonal-48 x 0.2ml

SP - COMPATIBILE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml

SP - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml

SP - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml