

Piastra fully-skirted 96x0.1ml (LP) - bianche

Cf 25pz - Piastra fully-skirted 96x0.1ml (LP) - bianche. Cf 25pz

**Product codes:**

Reference: 12 893 5681

Product gallery:**Product description:**

Piastra low profile da 0.1ml fully skirted. Utilizzabile su sistemi per PCR che utilizzano blocchi sia da 0.1ml che da 0.2ml, per PCR e PCR Real Time.

Questa piastra è robusta, impilabile e robotic friendly e non si incastrano quando vengono impilate.

- Colore bianco
 - Low Profile
 - Fully skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TIPO: piastre PCR
PSC - VOLUME: 0.1ml
PSC - FORMA: skirted
PSC - APERTURA: normale

PSC - POZZETTI: 96
PSC - FONDO: a V
PSC - CHIUSURA: tramite pellicola
PSC - CHIUSURA: tramite tappi
PSC - CHIUSURA: tappetino
PSC - APPLICAZIONI: PCR
PSC - APPLICAZIONI: biologia molecolare
PSC - APPLICAZIONI: PCR Real-Time FAST
PSC - APPLICAZIONI: PCR Real-Time
PSC - AUTOCLAVABILE: sì
PSC - COLORE: bianco
PSC - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITA': no
PSC - MATERIALE: PP
PSC - CONFORMITA': non soggetto a CE
PSC - IMPILABILE: sì
PSC - TEMPERATURA MAX: 100°C
PSC - TEMPERATURA MINIMA: -30°C
PSC - PACKAGE: bulk
PSC - PECULIARITA': piastra PCR 96 pozzetti da 0.1 per PCR e PCR Real Time e FAST skirted idonea per sistemi automatici - COLORE BIANCO
PSC - ACCESSORI OPZIONALI: film sigillanti
PSC - ACCESSORI OPZIONALI: tappetini sigillanti
PSC - ACCESSORI OPZIONALI: strip di tappi
PSC - IDONEO PER: stoccaggio
PSC - IDONEO PER: sistemi automatici
PSC - IDONEO PER: preparazione campioni
PSC - IDONEO PER: PCR Real-Time FAST
PSC - IDONEO PER: PCR
PSC - IDONEO PER: PCR Real-Time
PSC - PROFILO: basso tipo FAST/0.1
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro S 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Thermomixer C 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-100 96 x 0.2ml AgV
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-200 DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-220 DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler pro 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-221 DNA Engine Disciple Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-225 DNA Engine Tetrad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 96v Moto Alpha Unit 96 x 0.2ml
PSC - COMPATIBILE CON: MJ Research Inc. PTC-240 DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra T1 Thermocycler 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TGradient 96 x 0.2ml

PSC - COMPATIBILE CON: Whatman Biometra TProfessional 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TProfessional Standard Gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Whatman Biometra TRobot 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories C1000 Touch 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Connect Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories CFX96 Real-Time PCR Detection System 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 96v Alpha unit 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep realplex Silver 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient Silver 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler ep gradient 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96-77 96 x 0.2ml
PSC - COMPATIBILE CON: Eppendorf Mastercycler 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories T100 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Well Fast
PSC - COMPATIBILE CON: Bio-Rad Laboratories S1000 96-Deep Well
PSC - COMPATIBILE CON: Bio-Rad Laboratories PTC-100 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Tetrad 2 Chromo4 96 x 0.2ml
PSC - COMPATIBILE CON: Bio-Rad Laboratories DNA Engine Dyad Disciple 96v Alpha unit 96 x 0.2ml
PSC - CAPACITA': 0.1ml
PSC - PARETI: sottili 0.35mm
PSC - TRATTAMENTO SEGNALE: pareti bianche
PSC - FORZA SEGNALE: 180
PSC - POSSIBILE INTERFERENZA SEGNALE: inferiore a 2
PSC - RESISTENTE: congelamento a -20°C
PSC - EVAPORAZIONE CAMPIONE: minore del 3%
PSC - BINDING: DNA minore 0.3%
PSC - BINDING: proteine minore 0.8%
PSC - GRIGLIA ALFANUMERICA: si