

Strip 8x0.2ml (RP) extra-rigide - Ultra-trasparenti

Cf 120pz - Strip 8x0.2ml (RP) extra-rigide - Ultra-trasp. Cf 120pz

Product codes:

Reference: 12 954 3011

Product gallery:



Product description:

Strip da 8 microtubi compatibile con molti strumenti per PCR e PCR Real Time che accettano plastica per biologia molecolare a profilo regolare e alcuni modelli di sequenziatori

- Disegnata per applicazioni in PCR e PCR Real Time , compatibile anche con apparecchi come ABI / Life Technologies, Bio-Rad, Eppendorf e altri
- Striscia robusta non flessibile
- Versione profilo regolare da 0.2ml
- Colore neutro ultra trasparente
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

Product features:

SP - TIPO: strip da 8 microtubi

SP - TAPPO: senza

SP - VOLUME: 0.2ml

SP - CAPACITA': 0.2ml
SP - POZZETTI: 8
SP - PROFILO: regolare
SP - ALTEZZA: 20.9 / 21.8mm ca
SP - FORMA: striscia robusta
SP - FONDO: conico
SP - MATERIALE: PP
SP - APPLICAZIONI: PCR
SP - APPLICAZIONI: PCR Real Time
SP - AUTOCLAVABILE: si
SP - COLORE: neutro
SP - PUREZZA: DNA / RNA / DNase / RNase / Pyrogen FREE
SP - STERILITA': no
SP - CONFORMITA': non soggetto a CE
SP - TEMPERATURA MAX: 100°C
SP - TEMPERATURA MINIMA: -30°C
SP - PACKAGE: busta
SP - PECULIARITA': strip da 8x0.2 per PCR Real Time e PCR a striscia robusta - ultra trasparenti
SP - ACCESSORI OPZIONALI: griglia assemblaggio gialla
SP - IDONEO PER: PCR
SP - IDONEO PER: PCR Real Time
SP - PARETI: sottili 0.30mm
SP - TRATTAMENTO SEGNALE: nessuno
SP - FORZA SEGNALE: 100
SP - POSSIBILE INTERFERENZA SEGNALE: 100
SP - RESISTENTE: congelamento a -20°C
SP - EVAPORAZIONE CAMPIONE: minore del 3%
SP - BINDING: DNA minore 0.3%
SP - BINDING: proteine minore 0.8%
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000-96-Well Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Deep Well
SP - COMPATIBILE CON: Bio-Rad Laboratories-S1000-96-Well Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Thermomixer C-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV

SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TGradient-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TRobot-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-UNO II-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-310-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-3100-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-3130-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-3500 Genetic Analyzer-8-capillary
SP - COMPATIBILE CON: Applied Biosystems-3700-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-3730-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-ProFlex PCR System-32 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-ProFlex PCR System-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-SimpliAmp-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-Solid3-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-Veriti 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-Veriti Dx 96-Well Thermal Cycler-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml

SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
SP - COMPATIBILE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TOptical-96 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TPersonal-48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
SP - COMPATIBILE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml
SP - COMPATIBILE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
SP - COMPATIBILE CON: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-3500xL Genetic Analyzer-24-capillary
SP - COMPATIBILE CON: Applied Biosystems-ViiA 7 Real-Time PCR System-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-MiniAmp Plus Thermal Cycler-96 x 0.2ml
SP - COMPATIBILE CON: Applied Biosystems-MiniAmp Thermal Cycler-96 x 0.2ml