



# Blue Gel Loading Buffer with DNA Stain - 5 x 1,8ml

Blue Gel Loading Buffer with DNA Stain

Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain



**Product codes:**

Reference: JN PCR-274

**Product gallery:**

**Product description:**

Blue Gel Loading Buffer with DNA Stain - Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain

5 x 1,8ml

**Shipping Conditions:** shipped on gel packs

**Storage Conditions:** store at -20 °C  
store dark

**Shelf Life:** 12 months

**eClass Nummer:** 32160502

Product Web Link: [Product Web Link](#)

Product Data Sheet: [Download Data Sheet](#)

Product Safety Data Sheet: [Download Product Safety Data Sheet](#)

**Product features:**

RGT - TYP: Reagenti

RGT - TYP: Polymerase chain reaction PCR and reverse transcriptase polymerase chain reaction RT PCR products

RGT - ABLAUF: 12 months

RGT - TRANSPORT: shipped on gel packs



RGT - CAPACITA': 5 x 1

RGT - CONSERVAZIONE: Store at -20 °C

RGT - CAS NUMBER: -